

# AGED CARE WAREHOUSE CLEARANCE

ALL PRODUCTS  
BETWEEN  
**50-70% OFF**

HURRY,  
WHILST  
STOCK  
LASTS!

Item Number	Item Description	WAS	NOW	QTY
BEA007214	Alarm - Invisabeam - Invisalert - Pager	\$ 839.00	\$ 251.70	3
BEA007215	Alarm - Invisabeam - Remote with Power Supply	\$ 539.00	\$ 161.70	3
BEA009320	Bed Pole - Folding - Top Mounted - 32 x 50mm	\$ 249.00	\$ 74.70	2
BEA014605	Bed Cradle - Aspire - Fixed - Treated Steel ^	\$ 69.00	\$ 20.70	38
BEA024123	Bed End Set - CS7 wide - Amherst 42" w/ Ctrl Panel & T-burg - Solar Oak	\$ 869.00	\$ 260.70	1
BEA029814	Moulded ABS Bed Ends - Aidacare Ward Bed/FL250 ACUTE	\$ 260.00	\$ 78.00	1
BEA030050	Bed Mover/Transporter - Aidacare FL250 Bed (Set of 2))	\$ 809.00	\$ 404.50	1
BEA829110	Bed Side Rails - Full Length Beech - Aspire Homecare Bed	\$ 323.00	\$ 96.90	1
BEB045860	Bed - Aidacare AC3 V2 - King Single - 3 Function - Teak Bed Ends	\$ 2,950.00	\$ 1,475.00	2
BEM046725	Mattress - Aspire Lifecomfort Premier Care - Single - 180kg	\$ 869.00	\$ 434.50	5
BTC067108S	Shower Commode Elevating Legrest - Aspire - LEFT - Non-Tilt Commodes Only - Shiny Finish	\$ 370.00	\$ 111.00	1
BTC067110S	Shower Commode Elevating Legrest - Aspire - LEFT - Tilt in Space Commodes Only - Shiny Finish	\$ 370.00	\$ 111.00	1
BTC067115S	Shower Commode Elevating Legrest - Aspire - RIGHT - Tilt in Space Commodes Only - Shiny Finish	\$ 370.00	\$ 111.00	1
BTC069590	Shower Commode Seat - 4HC Verto TIS - Closed Front - 700mm Wide	\$ 739.00	\$ 221.70	3
BTC087630	Chest belt	\$ 109.00	\$ 32.70	1
BTT145110	Over Toilet Aid - Padded Toilet Seat only - Aspire - Grey	\$ 178.00	\$ 53.40	2
BTT146300	Over Toilet Aid - K Care - Wide 500mm 135Kg - Zinc	\$ 278.00	\$ 83.40	18
CHP198960	Chair Pressure Relief - Aspire COVE Chair - Grey/Andreas Pumice - Medium - 160kg ^	\$ 3,689.00	\$ 1,844.50	1
CHP198965	Chair Pressure Relief - Aspire COVE Chair - Navy/Bronx Steel - Medium - 160kg ^	\$ 3,689.00	\$ 1,844.50	1
CHP198970	Chair Pressure Relief - Aspire COVE Chair - Grey/Andreas Pumice - Large - 160kg ^	\$ 3,689.00	\$ 1,844.50	4
CHP198975	Chair Pressure Relief - Aspire COVE Chair - Navy/Bronx Steel - Large - 160kg	\$ 3,689.00	\$ 1,844.50	3
CHP212700	Hilite Chair Tray Adjustable - Grey/Slate	\$ 509.00	\$ 152.70	1
CHP216300	Chair Footrest - Princess Chair/Bed & Deluxe - Air Comfort	\$ 477.00	\$ 143.10	3
LAC353030	Linen Bag Waterproof Black	\$ 70.00	\$ 21.00	2
LIC354115	Soft Strap - P200 Weigh Scale	\$ 58.00	\$ 17.40	2
LIC355820	Limit Switch	\$ 31.30	\$ 9.39	2
LIC355840	Limit Switch Shaft	\$ 5.98	\$ 1.79	4
LIC355850	Limit Switch Roller	\$ 7.45	\$ 2.24	4
LIC359150	Carriage Hook for Multirall or Likorall	\$ 187.00	\$ 56.10	1
LIC359260	Extension Arm - Liko Multirall - Fits Carriage S65, Prod. No. 3136011.	\$ 194.00	\$ 58.20	2
LIC359370	Hand control 2 button - With display - LikoGuard	\$ 929.00	\$ 278.70	1
LIC359640	Carriage with Brake for Likorall R2R for 90deg turning of motor	\$ 270.00	\$ 81.00	1
LIC360670	HandControl HB80 - Liko	\$ 374.00	\$ 112.20	1
LIC360752	FreeSpan Rail 2m - with Likorall Trolley - 250kg	\$ 1,959.00	\$ 979.50	1
LSS382500	Lifter Sling - Full Body Polyester w/ Div Leg - Kerry - Large LR408N3	\$ 376.00	\$ 188.00	1
LSS390674	Integrated Weigh Scale - Allegro	\$ 2,679.00	\$ 803.70	18
LSS393360	RollOn - Max 160 kg (350 lbs).	\$ 2,800.00	\$ 1,400.00	1
LSS394160	Sling - SafetyVest Mod 93 - Liko - Med - w/4 Loop Clips Fits Sabina II 350 SlingBar	\$ 809.00	\$ 242.70	1
LSS394190	Sling - ComfortVest Mod 95 - Liko - Med - w/4 Loop Clips Fits Sabina II Comfort SlingBar.	\$ 799.00	\$ 239.70	1
LSS394200	Sling - ComfortVest Mod 95 - Liko - Lge - w/4 Loop Clips Fits Sabina II Comfort SlingBar.	\$ 829.00	\$ 248.70	1
LSS395190	Liftstraps XL (set) - Liko - Poly - Fits FlexoStretch and UltraStretch - 400kg/set SWL	\$ 639.00	\$ 191.70	1
MWS449610	WChair - Tilt in Space - Aspire Rehab RS - 460mm - 24in Wheels - 180kg - Gunmetal Grey	\$ 2,479.00	\$ 1,239.50	2
PTA494053	HipSaver - Slim Fit - Pant Only - Female - X Sml	\$ 65.00	\$ 45.50	3
PTA513800	HipSaver - SlimFit - High Compliance w/ Tailbone Pad - Female - Lge	\$ 170.00	\$ 119.00	1
PTA521913	Shoe - Sheepskin Hi Temp - DiabPro - Closed-Toe - Medium	\$ 210.00	\$ 147.00	2
PTA521918	Shoe - Sheepskin Hi Temp - Slimline Wrap-Around - Closed-Toe - Small	\$ 202.00	\$ 141.40	1
PTA534102L	Limb/Cast Protection - Limbo - Adult Below Elbow - Waterproof - Lge	\$ 61.00	\$ 42.70	1
PTC550300	Cushion - Evolight - Otto Bock - 55.5x48cm	\$ 659.00	\$ 329.50	1
PTC551290	Cushion - ROHO AirLite 14 x 16in	\$ 389.00	\$ 194.50	1
PTC571913	Cushion - JAY - Lite 16x18in	\$ 509.00	\$ 254.50	1
PTC573427	Cushion - JAY - Union 22 x 24in	\$ 779.00	\$ 389.50	1
PTC586900	Cushion - ROHO Nexus Spirit 16x18in	\$ 726.00	\$ 363.00	1
PTC590300	Cushion - ROHO Quadtro Select LP 11x11in (6x6cell)	\$ 1,119.00	\$ 559.50	1
PTC594921	Cushion - ROHO Hybrid Elite 1R 16x20in Sensor Ready	\$ 1,079.00	\$ 539.50	1
PTC599001	Cushion Solution - Pride Stealth - 18 x 16 in	\$ 579.00	\$ 289.50	6
PTC599201	Cushion Solution - Pride Stealth - 18 x 20 in	\$ 579.00	\$ 289.50	6
PTC608285	Cushion Cover - Domino - 3D - 16.5 x 18in	\$ 294.00	\$ 88.20	1
PTM638015	Mattress - Aspire HYBRID - King Single - 9 Cell - 1980x1050x170mm - 180kg (Pump Sep: PTM638030)	\$ 1,939.00	\$ 969.50	3
PTM639300T	Mattress Replacement - Aspire ACTIVE AIR 8 TURN - Single - 20 Cells - 2000x900x200mm - 200kg	\$ 2,821.00	\$ 1,410.50	1
PTS640022	JOBSTTravel Socks Black Size 2	\$ 58.00	\$ 17.40	1
PTS640150	Stockings - JOBST UltraSheer - Knee - 15-20mmHg - C/Toe - Med - Blk 1pair	\$ 70.00	\$ 21.00	3
PTS641975CLBEXL	Stockings - Venosan 6002 (Nyl/Lyc) - Knee - 23-32mmHg - Unisex - C/Toe - Long - XL - Bge	\$ 118.00	\$ 35.40	1
WAT900700	Medicine Cart - Aidacare - 650x480x980mm - White - 3 Drawers + 1 Dump Drawer	\$ 2,869.00	\$ 1,434.50	4
WAT905350-O	Isolation Cart Avalo 9H - DR003 w/ No Lock & RHS Handle - YELLOW (dividers incl)	\$ 2,374.00	\$ 1,187.00	1
WAT907550	Avalo Keys	\$ 25.10	\$ 7.53	1
WAT907665	Avalo Writing Surface - Light Cream	\$ 173.00	\$ 51.90	1
WAT909230	Cart Accessory - Aidacare/Tribute - Sharp Container	\$ 184.00	\$ 55.20	1