

WAREHOUSE CLEARANCE

ALL PRODUCTS BETWEEN **50-70% OFF**

HURRY, WHILST STOCK LASTS!

Item Number	Item Description	WAS	NOW	QTY
BEA829110	Bed Side Rails - Full Length Beech - Aspire Homecare Bed	\$ 323.00	\$ 96.90	92
BEA829150	Upholstered Surround Kit - Allure - Portobello Grey - Aidacare AC4 - Single	\$ 579.00	\$ 173.70	24
BEB045105	Bed - Jomor - Cuddle Bed Configuration - 4 Function - 200kg	\$ 13,969.00	\$ 4,190.70	1
BEB045150	Bed - UC - Floorline Single w/ Knee Bend	\$ 5,144.90	\$ 1,543.47	3
BEB045850	Bed - Aidacare AC3 V2 - Single - 3 Function - Teak Bed Ends	\$ 2,450.00	\$ 1,225.00	5
BEB050100	Hospital Bed - Malvestio VIVO ICU Bed w/Integrated Air Mattress - 2190-2420x1080x480-820mm - 260kg	\$ 43,879.00	\$ 21,939.50	1
BEB050201	Hospital Bed - Malvestio Botero Bariatric Ward Bed w/Integrated Weigh Scale - 400kg	\$ 21,989.00	\$ 10,994.50	1
BEB050310	Hospital Bed - Aspire Gynaecology Delivery Bed - Electric - 2175x1060x620-890mm - 175kg	\$ 9,495.00	\$ 4,747.50	1
BTC070900SP3	SWEP/DES - Shwr Commode SA Frest + (Paediatric) w/Closed Seat & Pan Carrier	\$ 2,519.00	\$ 755.70	4
BTS118200	Change Table - Mobile - Kerry - Electric	\$ 10,729.00	\$ 3,218.70	1
BTS152100	Swivel Bather - Freedom - HBA431 With Rotating Seat - 110kg	\$ 497.00	\$ 149.10	20
CHP197706	Chair Powerlift Recline - Aspire Signature - Space Saver - Small - Heritage - 145kg	\$ 2,099.00	\$ 1,049.50	17
CHP197741	Chair Powerlift Recline - Aspire Signature - Dual Action - Small - Heritage Vinyl - 155kg	\$ 2,269.00	\$ 1,134.50	19
CHP198965	Chair Pressure Relief - Aspire COVE Chair - Navy/Bronx Steel - Medium - 160kg ^	\$ 3,689.00	\$ 1,844.50	10
CHP198970	Chair Pressure Relief - Aspire COVE Chair - Grey/Andreas Pumice - Large - 160kg ^	\$ 3,689.00	\$ 1,844.50	3
CHP198975	Chair Pressure Relief - Aspire COVE Chair - Navy/Bronx Steel - Large - 160kg	\$ 3,689.00	\$ 1,844.50	4
CHP198976	Chair Pressure Relief - Aspire COVE Chair - Grey Vinyl - Medium - 160kg	\$ 3,519.00	\$ 1,759.50	18
CHP198976	Chair Pressure Relief - Aspire COVE Chair - Grey Vinyl - Medium - 160kg	\$ 3,519.00	\$ 1,759.50	18
LIC350055	PRISM TROLLEY RING GUIDE- 8 WHEEL TROLLEY - ABS (UK)	\$ 88.00	\$ 26.40	153
LIC350060	PRISM TROLLEY RING GUIDE - ABS PLASTIC (REHABASSIST)	\$ 57.00	\$ 17.10	447
WAF750150	Seat Walker - Freedom - Stride - Red - 135kg	\$ 369.00	\$ 184.50	105

Item Number	Item Description	WAS	NOW	QTY
BEA007444	Alarm - Cura - Dual Adaptor - Black	\$ 57.00	\$ 17.10	5
BEA007473	Alarm - Cura - Monitor Soft Plastic Mounting Holder (OLD - use with 2613, 2610, 2608)	\$ 40.00	\$ 12.00	1
BEA007496	Alarm - Cura - 30 Channel Display Pad	\$ 410.00	\$ 123.00	1
BEA007514	Alarm - 6mm Extension Tele Plug - To Suit Nurse Call Systems	\$ 30.70	\$ 9.21	2
BEA007520	Alarm - Monitor - Mono - To Suit Nurse Call Systems	\$ 184.00	\$ 55.20	5
BEA007522	Alarm - Cura - Cable Connector - 6.35 to 3.5 Mono (T&S) - 2m	\$ 31.50	\$ 9.45	4
BEA007523	Alarm - Cura - Cable Connector - 6.35 Stereo (T&S) - 2m	\$ 31.50	\$ 9.45	1
BEA007536	Alarm - Cura - Pad Connector - USB to RJ-12 Phone Plug (Wired Pin 2&5) - 2m	\$ 31.50	\$ 9.45	5
BEA007681	Alarm - Alerta - Mains Power Adaptor - Detect Receiver ONLY	\$ 48.60	\$ 14.58	5
BEA007682	Alarm - Alerta - Mains Power Adaptor - Detect Motion ONLY	\$ 48.60	\$ 14.58	5
BEA007850	Alarm - Sofihub TEQ Home (1 x SOFIHUB TEQ Home + 8 Sensors)	\$ 1,419.00	\$ 709.50	6
BEA014605	Bed Cradle - Aspire - Fixed - Treated Steel ^	\$ 69.00	\$ 20.70	17
BEA015000	Overbed/Chair Table - K Care - Tilt Dual HD - Wood/Brown	\$ 493.00	\$ 147.90	1
BEA016571	Headboard - iCare - Long single - Black Vinyl	\$ 579.00	\$ 173.70	1
BEA019950	Wonder Sheet - Fitted - Single Size - 45cm Depth - ECRU	\$ 124.00	\$ 62.00	6
BEA024120	Bed End Set - CS7 wide - Drake 42" w/ Ctrl Panel - Biltmore Cherry	\$ 425.30	\$ 127.59	2
BEA024121	Bed End Set - CS7 wide - Drake 42" w/ Ctrl Panel & T-burg - Biltmore Cherry	\$ 789.00	\$ 236.70	5
BEA024123	Bed End Set - CS7 wide - Amherst 42" w/ Ctrl Panel & T-burg - Solar Oak	\$ 869.00	\$ 260.70	2
BEA024124	Assist Bar - Invacare - for Invacare CS7 Wide Deck Bed	\$ 330.00	\$ 99.00	1
BEA024128	Assist Rail - Invacare - for Invacare CS7 Wide Deck Bed	\$ 402.00	\$ 120.60	2
BEA024131	Self Help Pole - Invacare - for Invacare All CS Beds	\$ 509.00	\$ 152.70	4
BEA028100	Overbed/Chair Table - Aspire - Wide C Base - Thermoform Recessed Split Tilt Top - Cream	\$ 349.00	\$ 174.50	2
BEA029813	Mounting Bracket ONLY (NO FL250 HARDWARE) Bed Ends - Aidacare FL250 ACUTE	\$ 312.00	\$ 93.60	2
BEA029814	Moulded ABS Bed Ends - Aidacare Ward Bed/FL250 ACUTE	\$ 260.00	\$ 78.00	14
BEA029845	Bed End Shrouds - Beech - Single - Aspire COMMUNITY CARE Bed (Premium)	\$ 300.00	\$ 90.00	4
BEA029866	Bed Side Rails - Extended Length Beech - Aspire PREMIUM HOMECARE Bed	\$ 366.00	\$ 109.80	1
BEA029871	Bed Transport Kit - Aspire LIFESTYLE Community Bed	\$ 103.00	\$ 30.90	31
BEA030050	Bed Mover/Transporter - Aidacare FL250 Bed (Set of 2))	\$ 809.00	\$ 404.50	3
BEA824130	Bed Moving Trolley - Generic	\$ 1,029.00	\$ 308.70	5
BEA826290	Bedside Cabinet - Hansel - Timber (choose stain colour)	\$ 879.00	\$ 439.50	2
BEA828422	Bedside Cabinet - Belmont Timber (RAW) 2Drws & O/Space - Classic Oak	\$ 819.00	\$ 409.50	7
BEA829110	Bed Side Rails - Full Length Beech - Aspire Homecare Bed	\$ 323.00	\$ 96.90	92
BEA829150	Upholstered Surround Kit - Allure - Portobello Grey - Aidacare AC4 - Single	\$ 579.00	\$ 173.70	24
BEB045100	Bed - Long Double - 4 Function - 400kg	\$ 8,719.00	\$ 2,615.70	1
BEB045105	Bed - Jomor - Cuddle Bed Configuration - 4 Function - 200kg	\$ 13,969.00	\$ 4,190.70	1
BEB045140	Bed - UC - Floorline Single	\$ 4,395.00	\$ 1,318.50	1
BEB045150	Bed - UC - Floorline Single w/ Knee Bend	\$ 5,144.90	\$ 1,543.47	3
BEB045567	Bed - Invacare - CS7-42" Wide Deck (bed ends sep)	\$ 5,659.00	\$ 2,829.50	1
BEB045569	Bed - Invacare - CS5 without head and foot board	\$ 3,709.00	\$ 1,854.50	1
BEB045850	Bed - Aidacare AC3 V2 - Single - 3 Function - Teak Bed Ends	\$ 2,450.00	\$ 1,225.00	5
BEB045860	Bed - Aidacare AC3 V2 - King Single - 3 Function - Teak Bed Ends	\$ 2,950.00	\$ 1,475.00	1
BEB045905	Bed - Aidacare AC3 King Single - 3 Function - CL - Beech Head & Foot Boards - 220kg	\$ 3,409.00	\$ 1,704.50	1
BEB046000	Bed - Aidacare Ward Bed - 4 Function - Single - CL - 250kg (Siderails Not Included)	\$ 3,395.00	\$ 1,697.50	1
BEB046295	Bed End Set Custom (Tallo - Timber (RAW) - With Control Panel Cut-out) - FL250/FLX - King Single	\$ 889.00	\$ 266.70	4
BEB046605	Reprogramming Handset - Safety Egress Height Adjustment - Aidacare FL250/FLX Bed	\$ 200.00	\$ 60.00	15
BEB050100	Hospital Bed - Malvestio VIVO ICU Bed w/Integrated Air Mattress - 2190-2420x1080x480-820mm - 260kg	\$ 43,879.00	\$ 21,939.50	1
BEB050201	Hospital Bed - Malvestio Botero Bariatric Ward Bed w/Integrated Weigh Scale - 400kg	\$ 21,989.00	\$ 10,994.50	1
BEB050310	Hospital Bed - Aspire Gynaecology Delivery Bed - Electric - 2175x1060x620-890mm - 175kg	\$ 9,495.00	\$ 4,747.50	1
BEC046071	Bedside Commode - K Care - Economy - 440mm Champagne	\$ 324.00	\$ 162.00	1

WAREHOUSE CLEARANCE

ALL PRODUCTS BETWEEN **50-70% OFF**

HURRY, WHILST STOCK LASTS!

BEC046130	Bedside Commode - K Care - Deluxe - 440mm Mulberry - Zinc	\$ 661.90	\$ 330.95	2
BEM046725	Mattress - Aspire Lifecomfort Premier Care - Single - 180kg	\$ 869.00	\$ 434.50	5
BEM046736	Mattress - Aspire Lifecomfort Premier Acute - Ward Bed Single - 250kg	\$ 989.00	\$ 494.50	9
BEM046810	Mattress - Lifecomfort 1930x850x130mm (Cover Only)	\$ 257.00	\$ 128.50	1
BEP054300	Pillow Cover - Pyramid Ruffled Slip only	\$ 54.00	\$ 27.00	4
BEP054710	Pillow - Cuddle Up	\$ 146.00	\$ 73.00	1
BEP054760	Pillow Cover - Lucky 7 Body Bed Acc - Pillow	\$ 58.00	\$ 29.00	1
BTC062500	Shower Commode Seat - K Care - H/Made - Closed 550mm	\$ 839.00	\$ 251.70	1
BTC062700	Shower Commode Seat - K Care - H/Made - Closed 650mm	\$ 1,029.00	\$ 308.70	3
BTC062900	Shower Commode Seat - K Care - H/Made - Open 550mm - Blue	\$ 579.00	\$ 173.70	1
BTC064100	Shower Commode Seat - Pressure-Care - K Care - H/Made - Open 600mm	\$ 1,059.00	\$ 317.70	1
BTC066300	Bedpan and Bowl Carrier - 4HC (& KCare) - for Shower commodes	\$ 210.00	\$ 63.00	2
BTC067103	IV Pole - Adjustable - Aspire Shower Commodes	\$ 276.00	\$ 82.80	1
BTC067108S	Shower Commode Elevating Legrest - Aspire - LEFT - Non-Tilt Commodes Only - Shiny Finish	\$ 370.00	\$ 111.00	2
BTC067109S	Shower Commode Elevating Legrest - Aspire - RIGHT - Non-Tilt Commodes Only - Shiny Finish	\$ 370.00	\$ 111.00	2
BTC069590	Shower Commode Seat - 4HC Verto TIS - Closed Front - 700mm Wide	\$ 739.00	\$ 221.70	1
BTC069625	Pan & Bowl Carrier - Suits 4HC 550mm width Commodes	\$ 1.00	\$ 0.30	3
BTC070600SP1	SWEP/DES - Shwr Commode Transit Low w/Open Seat & Pan Carrier	\$ 1,199.00	\$ 359.70	7
BTC070900SP3	SWEP/DES - Shwr Commode SA Frest + (Paediatric) w/Closed Seat & Pan Carrier	\$ 2,519.00	\$ 755.70	4
BTC071350	Shower Commode - K Care - Tilt - 500mm - S/Steel (seat sep)	\$ 6,649.00	\$ 1,994.70	1
BTC074000	Shower Commode Seat - K Care - H/Made - Open Left Side	\$ 899.00	\$ 269.70	1
BTC075300	Backrest - K Care - Curved Padded - STD Width - for K Care Shower Commodes - Fitted	\$ 310.80	\$ 93.24	1
BTC076600	Bowl Carrier - for K Care Shower Commodes - S/Steel	\$ 168.00	\$ 50.40	3
BTC077200	Pelvic/Chest Strap - for K Care Shower Commodes - Not Fitted - Requires Rivetting	\$ 81.70	\$ 24.51	2
BTC077900SP1	SWEP/DES - Shwr Commode Self Prop Low w/Open Seat & Pan Carrier	\$ 1,479.00	\$ 443.70	2
BTC077900SP2	SWEP/DES - Shwr Commode Self Prop Low w/Closed Seat & Pan Carrier	\$ 1,479.00	\$ 443.70	1
BTC079600	Contoured Headrest - for K Care Shower Commodes	\$ 999.00	\$ 299.70	1
BTC080360	Foot Board - K Care - Shower Commode - Grey Plastic	\$ 124.00	\$ 37.20	3
BTC085100	Bowl Carrier - K Care - 320mm - for K Care Shower Commodes - S/Steel	\$ 354.00	\$ 106.20	6
BTC087637	Headrest Bracket - for Aquatec Ocean VIP Shower Commodes	\$ 16.90	\$ 5.07	1
BTC087692	Comfort soft backrest Ocean Ergo XL	\$ 74.00	\$ 22.20	1
BTC089124	Shower Commode Seat - Aspire - Hand Made Pressure Reducing - 320mm - Closed Front	\$ 346.00	\$ 103.80	1
BTC089180	Shower Commode Seat - Aspire - Hand Made Pressure Reducing - 500mm - Closed Front	\$ 461.00	\$ 138.30	2
BTR101270	Grab Rail Economy 25mm Length 450mm S.S.	\$ 62.00	\$ 18.60	2
BTR101300	Grab Rail Economy 25mm Length 600mm S.S.	\$ 69.00	\$ 20.70	10
BTR101982	Grab Rail - Concealed Flange 32mm - 1200 x 32mm - White	\$ 104.00	\$ 31.20	1
BTR102122	Grab Rail Anti Slip (Peened) Stainless Steel 450mm x 32mm Exposed Fix.	\$ 69.00	\$ 20.70	8
BTR103140	Grab Rail - Shower Rail - Waist shower grab rail - Brushed S/Steel	\$ 249.00	\$ 74.70	1
BTR103170	Grabrail Stainless Steel - L Shape with 90 degree return - Right hand	\$ 257.00	\$ 77.10	1
BTS105400	Bath Board Alton Bariatric 190kg	\$ 285.00	\$ 85.50	2
BTS106000	Shower Stool w/ Arms & Padded Seat - K Care - 550mm Seat - Zinc	\$ 778.40	\$ 233.52	1
BTS106510	Bath Stool - Breezy Every Day - 600mm - 225kg	\$ 107.00	\$ 32.10	1
BTS106700	Bath Chair - Myco Swivel Bather - CORNER	\$ 927.00	\$ 278.10	1
BTS107400	Hand Held Shower - Switchcock - S/Steel	\$ 45.00	\$ 13.50	1
BTS107800	Shower Hose Push On Dual Outlet (S.M.S.)	\$ 31.00	\$ 9.30	7
BTS108305	Bath Transfer Bench - Extension Legs	\$ 67.00	\$ 20.10	1
BTS114900HD	Heavy Duty Shower Chair w/ Arms - K Care - 440mm Seat - Zinc	\$ 170.00	\$ 51.00	11
BTS115200HD	Heavy Duty Shower Stool w/ Arms - K Care - 350mm Seat - Zinc	\$ 170.00	\$ 51.00	8
BTS115310	Shower Stool w/ Arms - K Care - 425mm Plastic Seat - Zinc-100kg	\$ 99.00	\$ 29.70	22
BTS116100	Shower Chair - Wall Mounted - Fold Away - K Care - AS Large Padded 125kg	\$ 2,358.57	\$ 707.57	2
BTS117700	Shower Stool w/ Arms & Padded Back & Seat - 125kg	\$ 319.00	\$ 95.70	3
BTS118200	Change Table - Mobile - Kerry - Electric	\$ 10,729.00	\$ 3,218.70	1
BTS118410	Change Table Side Rail - Kerry - Electric - RIGHT	\$ 1,479.00	\$ 443.70	1
BTS118411	Change Table Side Rail - Kerry - Electric - LEFT	\$ 1,369.00	\$ 410.70	1
BTS119300	Shower Stool - in U Frame Adjustable	\$ 97.00	\$ 29.10	2
BTS147130	Shower trolley hose kit	\$ 245.00	\$ 73.50	1
BTS148095	Shower Trolley LINER, PILLOW and FITTINGS- Aspire	\$ 1,209.00	\$ 362.70	1
BTS150100	Shower Chair - Freedom - HBA407 (Oval Tube) - Aluminium HD Shower Chair - 540mm - 250kg	\$ 275.00	\$ 82.50	1
BTS151110	Shower Chair - Saljol - Spa - Medium - Dark Grey - 150kg	\$ 220.00	\$ 66.00	2
BTS151120	Shower Chair - Saljol - Spa - Small - Dark Grey - 150kg	\$ 220.00	\$ 66.00	3
BTS152100	Swivel Bather - Freedom - HBA431 With Rotating Seat - 110kg	\$ 497.00	\$ 149.10	20
BTT146300	Over Toilet Aid - K Care - Wide 500mm 135Kg - Zinc	\$ 278.00	\$ 83.40	2
BTT147450	Chemical Spill KIT - Vikan Bucket with Absorb	\$ 215.00	\$ 64.50	2
BTT148900	Urinal Non-Spill Spout	\$ 29.90	\$ 8.97	5
BTT153200	Over Toilet Aid - K Care - 445mm - Zinc	\$ 132.00	\$ 39.60	5
BTT153200HD	Heavy Duty Over Toilet Aid w/splash guard - 445mm - Zinc	\$ 161.00	\$ 48.30	2
BTT153201	Over Toilet Aid w/Splash Guard - CareQuip - 440mm - SWL 150kg	\$ 149.00	\$ 44.70	31
BTT157118	Toilet Grab Rail - Throne Spacer Bracket - S/Steel	\$ 59.00	\$ 17.70	3
BTT163110	Over Toilet Aid - Freedom - HBA370-HD - Aluminium HD Dignity Split Seat Chair - 560mm - 200kg	\$ 325.00	\$ 162.50	2
CHP192600	Chair Lowback - Westgate Bridge Chair - 450mm - pc steel	\$ 267.00	\$ 133.50	1
CHP192900	Chair Highback - Comfort Chair Deluxe - 580mm - pc steel	\$ 629.00	\$ 314.50	4
CHP193000	Chair Footrest - Adjustable - Comfort/Harbour/Westgate Chairs	\$ 184.00	\$ 92.00	6
CHP194550	Lounge Chair - Cornwall Std. Fabric or Vinyl	\$ 1,229.00	\$ 614.50	1
CHP197706	Chair Powerlift Recline - Aspire Signature - Space Saver - Small - Heritage - 145kg	\$ 2,099.00	\$ 1,049.50	17
CHP197715	Chair Powerlift Recline - Aspire Signature - Space Saver - Medium - Lagoon - 145kg	\$ 2,099.00	\$ 1,049.50	1
CHP197716	Chair Powerlift Recline - Aspire Signature - Space Saver - Medium - Heritage Vinyl - 145kg	\$ 2,099.00	\$ 1,049.50	2

WAREHOUSE CLEARANCE

ALL PRODUCTS BETWEEN **50-70% OFF**

HURRY, WHILST STOCK LASTS!

CHP197725	Chair Powerlift Recline - Aspire Signature - Space Saver - Large - Lagoon - 145kg	\$ 2,099.00	\$ 1,049.50	2
CHP197735	Chair Powerlift Recline - Aspire Signature - Dual Action - Small - Mink - 155kg	\$ 2,269.00	\$ 1,134.50	1
CHP197740	Chair Powerlift Recline - Aspire Signature - Dual Action - Small - Lagoon - 155kg	\$ 2,269.00	\$ 1,134.50	1
CHP197741	Chair Powerlift Recline - Aspire Signature - Dual Action - Small - Heritage Vinyl - 155kg	\$ 2,269.00	\$ 1,134.50	19
CHP197770	Head/Armrest Cover Set - Aspire Signature Recliners - Mink	\$ 46.30	\$ 13.89	1
CHP197825	Variable Angle Leg Rest - Atlanta & Monaco Chairs	\$ 519.00	\$ 155.70	1
CHP197840	Armrest Cover w/ 2in Padding - Phoenix & Sorrento Chairs	\$ 659.00	\$ 197.70	1
CHP197862	Armrest Cover w/ 2in Padding - Monaco Chair	\$ 659.00	\$ 197.70	1
CHP197865	Leg Laterals - Phoenix & Sorrento Chairs	\$ 639.00	\$ 191.70	1
CHP197895	Removable Seat Cover - All Seating Matters Chairs	\$ 218.00	\$ 65.40	1
CHP198086	Chair Accessory - Configura - R/R Depth Adjustment Kit	\$ 314.00	\$ 94.20	1
CHP198088	Chair Accessory - Configura - Lateral Support Backrest Medium	\$ 449.00	\$ 134.70	1
CHP198810	Overnight Chair/Bed - Aspire SLEEPER Chair - Beige Vinyl - 150kg	\$ 3,109.00	\$ 1,554.50	3
CHP198820	Overnight Recliner Chair - Aspire RECOVERY Chair - Beige Vinyl - 150kg	\$ 1,109.00	\$ 554.50	1
CHP198958	Shoalhaven Tray Brackets	\$ 184.00	\$ 55.20	1
CHP198960	Chair Pressure Relief - Aspire COVE Chair - Grey/Andreas Pumice - Medium - 160kg ^	\$ 3,689.00	\$ 1,844.50	2
CHP198965	Chair Pressure Relief - Aspire COVE Chair - Navy/Bronx Steel - Medium - 160kg ^	\$ 3,689.00	\$ 1,844.50	10
CHP198970	Chair Pressure Relief - Aspire COVE Chair - Grey/Andreas Pumice - Large - 160kg ^	\$ 3,689.00	\$ 1,844.50	3
CHP198975	Chair Pressure Relief - Aspire COVE Chair - Navy/Bronx Steel - Large - 160kg	\$ 3,689.00	\$ 1,844.50	4
CHP198976	Chair Pressure Relief - Aspire COVE Chair - Grey Vinyl - Medium - 160kg	\$ 3,519.00	\$ 1,759.50	18
CHP198991	Advanced Immersion Backrest Cushion - Medium - Navy - Aspire COVE Chair	\$ 427.00	\$ 128.10	4
CHP198992	Advanced Immersion Backrest Cushion - Large - Grey - Aspire COVE Chair	\$ 427.00	\$ 128.10	2
CHP204500	Chair Highback - K Care Hilite - 490mm - Grey/Slate	\$ 560.00	\$ 168.00	3
CHP205000	Hilite Chair Mulberry	\$ 560.00	\$ 168.00	1
CHP208500	Chair Lowback - K Care Utility - 440mm - Vanilla/Champ	\$ 355.00	\$ 106.50	7
CHP211100	Hilite Chair Mobile Wheels Vanilla	\$ 709.00	\$ 212.70	2
CHP211400	Hilite Chair Neckrest - Grey/Slate	\$ 436.00	\$ 130.80	1
CHP216300	Chair Footrest - Princess Chair/Bed & Deluxe - Air Comfort	\$ 477.00	\$ 143.10	3
CHP216352	Chair Footrest - Air Comfort - Deluxe V2 - Large old style	\$ 331.00	\$ 99.30	3
CHP216353	Chair Tray - Air Comfort - Deluxe V2 - Large	\$ 385.00	\$ 115.50	1
CHP216400	Chair Tray Small - Princess Chair/Bed - Air Comfort	\$ 149.00	\$ 44.70	4
CHP217200	Chair Footrest - Air Comfort - My Comfort Chair	\$ 599.00	\$ 179.70	1
CHP227100	Chair Powerlift Recline - Aspire Verona - Mushroom - Single Motor - 130kg ^	\$ 1,659.00	\$ 829.50	5
CHP227115	Pressure Reducing Overlay - Aspire Maui/Verona Lift Recline Chairs	\$ 249.00	\$ 74.70	2
CHP227715	Oncology Armrest - Aspire ALTITUDE Single Action Chair	\$ 346.00	\$ 103.80	6
CHT100100	Treatment Chair - Greiner Multiline Next IT - Navy Blue Vinyl - 220kg	\$ 6,500.00	\$ 1,950.00	1
CHT100110	Treatment Chair - Greiner Multiline Next IT V2 - Navy Blue Vinyl - 220kg	\$ 8,249.00	\$ 2,474.70	1
CHT100200	Treatment Chair - Greiner Multiline Next DC - Light Grey Vinyl - 220kg	\$ 8,789.00	\$ 2,636.70	1
CHT100210	Treatment Chair - Greiner Multiline Next DC - Navy Blue Vinyl - 220kg	\$ 8,789.00	\$ 2,636.70	1
CHT100400	Treatment Chair - Aspire DIALYSIS Chair - Beige Vinyl - 150kg	\$ 1,869.00	\$ 560.70	2
CON228100	Bed Pad ComfortShield Deluxe Single Fitted 92x188cm	\$ 68.00	\$ 20.40	1
CON230400	Chair Pad - Buddies - Rectangle 600x500mm - Floral - 1500ml	\$ 39.00	\$ 11.70	1
CON240300	Cello AMD slip cotton feel Medium (Normal Abs) (Pkt 20) (1900ml)	\$ 48.60	\$ 14.58	4
CON241100	Cello EcoMidi (Pkt 20)	\$ 16.20	\$ 4.86	5
CON242400	TENA Slip Plus Small (Pkt 30)	\$ 59.70	\$ 17.91	3
CON243400	TENA Pad Lady Extra (Pkt 20)	\$ 34.40	\$ 10.32	3
CON244250	TENA Active Long (Pkt 24)	\$ 15.00	\$ 4.50	2
CON245714	Undergarment - Conni - Chantilly Ladies - Sizes 10-20 (Beige only)	\$ 41.02	\$ 12.31	2
CON245716	Undergarment - Conni - Classic Ladies - Sizes 10-20 (Beige only)	\$ 35.67	\$ 10.70	1
CON245757	Briefs Ladies - Conni Chantilly - Size16 (US 12) - Beige	\$ 41.20	\$ 12.36	2
CON248000	TENA Super Shields (Pkt 30)	\$ 22.40	\$ 6.72	1
CON250002	TENA Pants for Men Level 4 (Pkt 8)	\$ 34.60	\$ 10.38	11
CON250100	TENA Pants Mesh Fix Small	\$ 62.00	\$ 18.60	3
CON250200	Tena Fix Stretch Pants Medium (pk25)	\$ 71.00	\$ 21.30	2
CON250400	TENA Pants Mesh Fix Large (Pkt 25)	\$ 92.00	\$ 27.60	1
CON253000	MOLICARE Prem Men Pad 2D(Pkt 14)	\$ 19.40	\$ 5.82	8
CON257900	Briefs All In One Large 105-119cm (S.M.S.)	\$ 42.00	\$ 12.60	4
DCO100225	Dementia Care - Home Safety - Cupboard C-Lock	\$ 24.90	\$ 7.47	1
DCS100135	Dementia Care - Orientation Signage Bedroom Foam Board, Blue.	\$ 83.00	\$ 24.90	8
DCS100180	Dementia Care - Orientation Signs - Kitchen - Foam Board, Red	\$ 71.30	\$ 21.39	3
DCS100240	Dementia Care - Orientation Signs - Toilet - Foam Board, Yellow	\$ 71.30	\$ 21.39	1
DCS100270	Dementia Care - Orientation Signs - Wardrobe - Medium	\$ 47.40	\$ 14.22	3
DIN265490	Food Tray 550x400mm (S.M.S.)	\$ 31.00	\$ 9.30	4
DIN265740	Jug with Lid 1 Litre Natural	\$ 13.30	\$ 3.99	9
DIN265830	Tumbler Autoplas 230ml Natural	\$ 5.49	\$ 1.65	2
DIN266600	Tumbler Autoplas 230ml Green	\$ 5.49	\$ 1.65	6
DLD269200	Shoelace - Soesi Laces Black - 450mm 3 pairs	\$ 13.00	\$ 3.90	1
DLD269330	Comb - Long Handled Fine Comb	\$ 42.50	\$ 12.75	1
DLD269395	Leg Lifter - Double Handed	\$ 32.70	\$ 9.81	2
DLG271070	Card Holder (S.M.S.)	\$ 13.00	\$ 3.90	1
DLG272215	Lighting - Outdoor Sensor Light	\$ 71.00	\$ 21.30	7
DLK275030	Measuring 2 Cup Good Grips Angled	\$ 16.30	\$ 4.89	6
DLK275125	Smooth Potato Masher Good Grips	\$ 39.60	\$ 11.88	10
DLK275240	Cutlery - Rocker - Knife - Knife with Solid Handle	\$ 49.00	\$ 14.70	2
DLK275630	Peeler - Gordon with Clamp	\$ 194.00	\$ 58.20	1
DLK276380	Cutlery - Amefa - Contoured Knife Angled	\$ 83.69	\$ 25.11	1

WAREHOUSE CLEARANCE

ALL PRODUCTS BETWEEN **50-70% OFF**

HURRY, WHILST STOCK LASTS!

DLK276390	Cutlery - Amefa - Contoured Fork Angled Left	\$ 52.31	\$ 15.69	1
DLK276410	Cutlery - Amefa - Contoured Spoon Angled Left	\$ 93.00	\$ 27.90	1
DLK276580	Cups & Mugs - Feeding Cup & Adjustable Spout 8mm	\$ 34.10	\$ 10.23	1
DLK278340	Neater Eater - Soup Bowl & Saucer	\$ 95.00	\$ 28.50	1
FBS280195	Wrist Splint - Rolyan - D-ring Thermal - Spica Splint - Left - Med	\$ 90.00	\$ 27.00	1
FBS280260	Knee Brace - Immobiliser - 3 Panel 50cm	\$ 66.00	\$ 19.80	2
FBS280904	Ankle Brace - McDavid - Lace Up Ultralight - Blk - Lge	\$ 100.00	\$ 30.00	1
FBS281020	Shin Splint Wrap - Thermal - Donjoy - One Size	\$ 66.00	\$ 19.80	1
FBS281720	Soft Collar Small	\$ 16.20	\$ 4.86	6
FBS283350	Hand - Gloves - Take Off - Half Finger Large	\$ 121.00	\$ 36.30	2
FBS283630	Back - BodyAssist - Sacro Cynch Belt - Bge - Sml	\$ 113.00	\$ 33.90	1
FBS283700	Abdominal Support - Belt - BodyAssist - Large	\$ 73.00	\$ 21.90	2
FBS283810	Wrist/Thumb Brace - BodyAssist - Thermal - LGE/XLG	\$ 46.20	\$ 13.86	1
FBS283830	Squiz Ball NH Medium Rating	\$ 19.90	\$ 5.97	3
FBS284370	Knee Brace - Body Assist - Elastic Slip-on - XL	\$ 25.00	\$ 7.50	1
FBS284380	Thigh Support - BodyAssist - Elastic Slip on - Sml	\$ 57.00	\$ 17.10	2
FBS285750	Cervical Collar Child 5x40cm	\$ 18.70	\$ 5.61	4
FBS285780	Abdominal Binder 9" Regular	\$ 64.00	\$ 19.20	1
FBS301214	Neck Support - Aspen Collar - Cervical Collar - Adult XTall	\$ 101.00	\$ 30.30	1
FBS301221	Neck Support - Miami-J Collar - Cervical Collar - one Collar + pads - Sml	\$ 235.00	\$ 70.50	1
FBS301464	Finger - Rolyan - Finger Extension Splint - Sml	\$ 78.00	\$ 23.40	1
FBS301464L	Finger - Rolyan - Finger Extension Splint - Lge	\$ 78.00	\$ 23.40	1
FBS301527	Knee Hinged - Rehband - Elastic Wrap Genu Direxa - O/Patella - XX Sml	\$ 252.00	\$ 75.60	1
FBS3015295	Knee Hinged - McDavid - PSII - Thermal Slip-on PSII - O/Patella - XXL	\$ 188.00	\$ 56.40	1
FBS301541	Knee Immobiliser - Post-op Knee Immobiliser - Med	\$ 50.00	\$ 15.00	20
FBS301572XL	Ankle Brace - Oppo - Lace-up w/ Straps & Stays - Total stability - XL	\$ 68.00	\$ 20.40	1
FBS30157697	Shoe - Cosmac - Post Op - Standard - Male - Lge	\$ 48.80	\$ 14.64	2
FBS301579	Shoe - Post Op - Darco - Med Surg Shoe - Women - Sml	\$ 51.00	\$ 15.30	1
FBS301606	Foot Support - Achilles Heel Protection Sleeve (UNIVERSAL)	\$ 69.00	\$ 20.70	1
FBS3016127	Hernia Aid - Scudotex - Hernia Support Belt w/ Pads - XXXL	\$ 105.00	\$ 31.50	1
FBS301840	Hand Orthoses - RMI BendEase - Left	\$ 121.00	\$ 36.30	1
FBS301853	Neck Collar - NeuroFlex RestorativeTM Kentucky Kollar - Small	\$ 245.00	\$ 73.50	2
FBS302010	Elbow Orthoses - NeuroFlex Restorative - Med.	\$ 276.00	\$ 82.80	2
FBS302012	Elbow Orthoses - NeuroFlex Restorative - Elbow Separator	\$ 112.00	\$ 33.60	7
FBS3020492	Foot Orthoses - RM Safe Boot II - Large	\$ 102.00	\$ 30.60	1
FBS302056	Neck Orthoses - NeuroFlex Restorative Safe Spine TLSO - Average/Regular	\$ 333.00	\$ 99.90	3
FBS306050	Donjoy 82-9619 Ergoform Thumb Support - Size 1	\$ 29.90	\$ 8.97	1
FBS307540	MEDI R.210 Thumb Support - Right - SML	\$ 66.00	\$ 19.80	1
FPF290384	Exercise - Weights - Theraband Soft Weights Green 2.0kg	\$ 55.00	\$ 16.50	2
FPF290460	Exercise - Weights - Hand Grips Extra Strong	\$ 20.00	\$ 6.00	1
PPM100190	Exercise Therapy - Hasomed - Power Supply Unit - for RehaMove2	\$ 145.00	\$ 43.50	1
FPP292166	Walker Boot - Air Walker Boot - Range of Motion Small	\$ 203.00	\$ 60.90	1
FPP299100	Exercise System Digi Flex Black 4.1kg	\$ 92.00	\$ 27.60	1
FPP303000	Exercise Ball - 45cm	\$ 55.00	\$ 16.50	2
FPP303300	Exercise Ball - 55cm	\$ 65.00	\$ 19.50	3
FPP303600	Exercise Ball - 75cm	\$ 79.00	\$ 23.70	2
FPP303800	Physio Balls 100cm	\$ 105.00	\$ 31.50	1
FPP309915	Exercise - Theraband - Tubing with Soft Grip Handles 121cm Silver Super Heavy	\$ 54.00	\$ 16.20	1
FTW109940	Footwear - Dr Comfort Crew Sock Black - Crew Sock Black Small 6011	\$ 29.90	\$ 8.97	4
GMD300310	Sphygmomanometer - Aneroid Erka Switch	\$ 166.00	\$ 49.80	1
GMD300355	Sphygmomanometer - Mercury Free	\$ 799.00	\$ 239.70	8
GMD332330	Sphygmomanometer Basic Aneroid Velcro Cuff	\$ 38.50	\$ 11.55	1
GMD332345	Otoscope Black	\$ 167.00	\$ 50.10	1
GMD332720	Sphygmomanometer Mercurial Floor (S.M.S.)	\$ 450.00	\$ 135.00	1
GMD332770	Sphygmomanometer Adult Cuff (S.M.S.)	\$ 10.00	\$ 3.00	2
GMD332780	Sphygmomanometer Obese / Thigh Cuff Velcro (S.M.S.)	\$ 32.00	\$ 9.60	2
GMD332930	Sphygmomanometer - Velcro Cuff only Adult Large	\$ 37.60	\$ 11.28	2
GMD332974	One Piece Cuff & Bladder 2 Tube Adult	\$ 32.20	\$ 9.66	2
GMD332975	Bladder Latex Free 2 Tube Obese	\$ 71.00	\$ 21.30	1
GMS336326	KTR - Single Button Silicone Pendant Mono Plug - Jack Plug (3m Lead)	\$ 121.00	\$ 36.30	2
GMS340951	Donning Aid - Easy Slide for O/Toe Stockings - Sml	\$ 48.00	\$ 14.40	1
GSS340180	Donut Cushion Gel 140x60x35mm	\$ 136.00	\$ 40.80	5
LAC353109	Linen Bag - Half Size (OHS) - Yellow	\$ 18.70	\$ 5.61	8
LAC353171	Washnet Standard 40 x 60 with Zip & Tag	\$ 15.30	\$ 4.59	1
LAC353780	Trolley - Kerry - Soiled Linen Trolley Single S/S w/ Foot Operated Lid	\$ 769.10	\$ 230.73	2
LAC353800	Trolley - Kerry - Soiled Linen Trolley Double No Bag & Epoxy Coated	\$ 669.00	\$ 200.70	1
LAC354050	Trolley - Kerry - Moist Linen Trolley Large Base Only, Spring-Loaded & Galvanized	\$ 719.00	\$ 215.70	1
LAC354060	Trolley - Kerry - Moist Linen Trolley Spring Base (option) - Large	\$ 728.40	\$ 218.52	1
LAC354080	Soiled Linen Trolley Bag Support Only Raised	\$ 168.00	\$ 50.40	7
LAC500555	Cleaning - Enviro Dusting Floor Mop, 350Mm,Oates	\$ 84.00	\$ 25.20	1
LAC500635	Soiled Linen Trolley Lid - Aidacare - Plastic w/ Mounting Hardware	\$ 43.70	\$ 13.11	2
LAC500850	Linen Skip - Mild steel Epoxy powder coated + Plastic twin wheel - triple	\$ 619.00	\$ 185.70	3
LIC350055	PRISM TROLLEY RING GUIDE- 8 WHEEL TROLLEY - ABS (UK)	\$ 88.00	\$ 26.40	153
LIC350060	PRISM TROLLEY RING GUIDE - ABS PLASTIC (REHABASSIST)	\$ 57.00	\$ 17.10	447
LIC353360	Ceiling Hoist - Prism - Bracket End Stop - Long w/Charger	\$ 416.00	\$ 124.80	1
LIC353535	Ceiling Hoist - Prism - Trolley for Portable Hoists - Guildman Track - No Swivel	\$ 370.00	\$ 111.00	1

WAREHOUSE CLEARANCE

ALL PRODUCTS BETWEEN **50-70% OFF**

HURRY, WHILST STOCK LASTS!

LIC353550	Ceiling Hoist - Prism - Trolley for Portable Hoists - Fixed Eyelet	\$ 101.00	\$ 30.30	9
LIC354250	Ceiling Hoist - Prism - 6 Hook Carry Bar	\$ 280.00	\$ 84.00	1
LIC355470	Charging Block Assy Xtra	\$ 34.90	\$ 10.47	1
LIC355480	Mascot Charger UK 9940- Prism - For Transactive XTRA	\$ 312.00	\$ 93.60	1
LIC355590	Middle Cover Insert	\$ 75.00	\$ 22.50	1
LIC355630	Led Holding Kit. Led Surround	\$ 16.30	\$ 4.89	1
LIC355640	Bi-Colour Led	\$ 24.60	\$ 7.38	1
LIC355815	Battery Connector Assy - Prism - P200 Ceiling Hoist	\$ 45.50	\$ 13.65	1
LIC355820	Limit Switch	\$ 31.30	\$ 9.39	1
LIC355871	Back Cover Only - Prism - For P200	\$ 37.20	\$ 11.16	12
LIC355873	Shut Off Switch - Prism - For P200	\$ 76.00	\$ 22.80	3
LIC355882	P200 Spring Arm	\$ 6.24	\$ 1.87	4
LIC355965	CP200 CP136 HANDSET DOCK	\$ 66.00	\$ 19.80	2
LIC356005	CP200 EXTERNAL MEMBRANE SWITCH LABEL (Reval)	\$ 182.00	\$ 54.60	1
LIC356025	CP200 ON-OFF SWITCH LOOM	\$ 63.00	\$ 18.90	1
LIC356030	CP200 Motor Assembly	\$ 320.00	\$ 96.00	3
LIC357540	Carry Bar QRS X	\$ 1,449.00	\$ 434.70	1
LIC359150	Carriage Hook for Multirall or Likorall	\$ 187.00	\$ 56.10	12
LIC359200	Quick-release Hook Bolt - Multirall	\$ 219.00	\$ 65.70	1
LIC359270	Extension Arm Adapter - Liko Multirall or Likorall - Fits adapter MR/LR	\$ 91.00	\$ 27.30	1
LIC359300	Ceiling Hoist - LikoGuard L - 272kg (600lbs) SWL (no accessories)	\$ 4,919.00	\$ 1,475.70	1
LIC359330	Spreader Bar - Liko SlingGuard 670 twin - Max 363kg (800lbs)	\$ 839.00	\$ 251.70	2
LIC359360	Hand control 2 button - Without display - LikoGuard	\$ 709.00	\$ 212.70	7
LIC359370	Hand control 2 button - With display - LikoGuard	\$ 929.00	\$ 278.70	3
LIC359380	Standard Carriage - LikoGuard	\$ 219.00	\$ 65.70	12
LIC359390	Wall Charger - LikoGuard, AU/NZ	\$ 391.00	\$ 117.30	5
LIC359410	In Rail Charge Carriage - LikoGuard - (with Step up card and cable) - Max 363kg (800lbs)	\$ 559.00	\$ 167.70	3
LIC359420	In-Rail-Charge Charger - LikoGuard, AU	\$ 462.00	\$ 138.60	2
LIC359560	Hand Control - Likorall ES-4MT (for Lifting and Transfer Motor)	\$ 483.00	\$ 144.90	1
LIC359570	Hand Control - Likorall ES-6MTS (for Lifting, Transfer Motor and Switches)	\$ 579.00	\$ 173.70	2
LIC359800	Battery Charger - Likorall 200 (24V/0,22 A)	\$ 245.00	\$ 73.50	1
LIC359884	2.2m Straight Rail H70 - White - Liko	\$ 287.00	\$ 86.10	1
LIC360100	Rail Cover, 5 pcs a 2.0m (79 in) - Liko	\$ 89.00	\$ 26.70	3
LIC360630	Side Rail Switch, electrical - White - Liko	\$ 1,749.00	\$ 524.70	2
LIC360704	Traverse Rail Carrier, Wide/Raised 100mm - Liko	\$ 809.00	\$ 242.70	1
LIC360770	FreeSpan Rail 3m - with Multirall Trolley - 250kg	\$ 2,279.00	\$ 1,139.50	2
LIC360900	UltraTwistTM Wide - incl, hand control. Max 500 kg (1100 lbs).	\$ 3,829.00	\$ 1,148.70	2
LSS363150	ArjoHuntleigh 4 PT SLING, PADDED,LGE	\$ 899.00	\$ 269.70	1
LSS367750	Sling Liko Original H/B Thin Net Poly Reinf Leg Support w Belt	\$ 442.00	\$ 221.00	1
LSS370100	Leg Harness Liko Medium/Large	\$ 224.00	\$ 112.00	1
LSS371615	Lifter Sling - Oxford - Quickfit Poly Padded - Large	\$ 459.00	\$ 229.50	1
LSS371698	Lifter Sling - Oxford - Quickfit Deluxe Poly - Small	\$ 217.00	\$ 108.50	1
LSS371701	Lifter Sling - Oxford - Quickfit Deluxe Padded Polyester - Lge	\$ 539.00	\$ 269.50	2
LSS373315	Sling - Standup - Roze - Velcro - Keyhole - Universal Size	\$ 327.00	\$ 163.50	1
LSS374305	Sling - General Purpose - Fabric - XL	\$ 327.00	\$ 163.50	1
LSS374909	Sling - General Purp Head Support - Fab - Small	\$ 327.00	\$ 163.50	1
LSS374925	Sling - General Purpose Paediatric - Medium	\$ 261.00	\$ 130.50	1
LSS376300	Sling - Supreme Hygiene - Fabric Large	\$ 327.00	\$ 163.50	4
LSS376400	Sling - Standup - Roze - Buckle - Keyhole - Large	\$ 261.00	\$ 130.50	1
LSS378900	Lifter Sling - Standing Lifter - Kerry - w/ Buckle Large	\$ 297.00	\$ 148.50	1
LSS379000	Lifter Sling - Seating Sling - Kerry - Large	\$ 342.00	\$ 171.00	1
LSS380700	Lifter Sling - Gen Purp Mesh - Kerry - Large LR406NE3	\$ 297.00	\$ 148.50	2
LSS381000	Lifter Sling - Gen Purp Polyester - Kerry - Medium	\$ 609.00	\$ 304.50	1
LSS382400	Lifter Sling - Full Body Polyester w/ Div Leg - Kerry - Medium LR408N2	\$ 376.00	\$ 188.00	3
LSS382500	Lifter Sling - Full Body Polyester w/ Div Leg - Kerry - Large LR408N3	\$ 376.00	\$ 188.00	1
LSS383200	Lifter Sling - Hammock Mesh - Kerry - Small LR409NE1	\$ 499.00	\$ 249.50	1
LSS386203	Lifter - Invacare - Adaption Kit - Jasmine	\$ 326.00	\$ 97.80	1
LSS388800	Limb Lifter - Single Patient Use - 850mm x 140mm - 120kg per Limb (Box of 10)	\$ 209.00	\$ 62.70	1
LSS388805	Limb Lifter - Single Patient Use - Wide - 850mm x 250mm - 120kg per Limb (Box of 10)	\$ 247.00	\$ 74.10	1
LSS388810	Pannus Support - SallySling - Single Patient Use - 120kg (Box of 5)	\$ 191.00	\$ 57.30	2
LSS390674	Integrated Weigh Scale - Allegro	\$ 2,679.00	\$ 803.70	2
LSS390917	Lifter Sling - Aspire - MED - Deluxe General Purpose w/Head Supp - CLIPS - Mesh - 200kg	\$ 423.00	\$ 211.50	3
LSS390945	Lifter Sling - Aspire - MED - Deluxe Hygiene/Access w/Head Supp - CLIPS - Poly - 200kg	\$ 480.00	\$ 240.00	1
LSS390965	Lifter Sling - Aspire - MED - All Day General Purpose w/Head Supp - CLIPS - Parasilk - 200kg	\$ 539.00	\$ 269.50	1
LSS393360	RollOn - Max 160 kg (350 lbs).	\$ 2,800.00	\$ 1,400.00	1
LSS393410	Wall Mounted Battery Charger, AU/NZ for CBL Battery.	\$ 679.00	\$ 203.70	1
LSS393480	Spreader Bar - Liko Universal SlingBar 350 - Max 300 kg (660 lbs).	\$ 353.00	\$ 105.90	1
LSS393540	Spreader Bar - Liko Universal SlingBar 600 - Max 300 kg (660 lbs).	\$ 490.00	\$ 147.00	1
LSS393700	Liko Adapter 12 mm - Fits Uno with serial no. from 31825, LikoLight & Viking - Max 300kg	\$ 44.90	\$ 13.47	3
LSS393800	Adapter 22 for Quick-release Hook - Liko - Fits Liko MasterLift from -89 and Likorall.	\$ 53.00	\$ 15.90	1
LSS394400	Sling - Liko Ultrasling Mod 6- Poly w/ reinforced leg suppts - short - 500kg	\$ 926.00	\$ 277.80	2
LSS394440	Sling - Liko SilhouetteSling Mod 22 - Mesh - Lge - 200kg	\$ 539.00	\$ 269.50	1
LSS394450	Sling - Liko SilhouetteSling Mod 22 - Mesh - XS - 200kg	\$ 435.00	\$ 217.50	3
LSS394460	Sling - Liko SilhouetteSling Mod 22 - Mesh - Sml - 200kg	\$ 453.00	\$ 226.50	4
LSS394500	Sling - Liko HygieneSling w/ Belt - Poly - Reinforced leg supports - Lge - 200kg	\$ 519.00	\$ 259.50	1
LSS394660	Sling - Liko Original H/B MOD 200 - Net Polyester - Reinforced Leg Supports - w/belt - XS - Child	\$ 519.00	\$ 259.50	1

WAREHOUSE CLEARANCE

ALL PRODUCTS BETWEEN **50-70% OFF**

HURRY, WHILST STOCK LASTS!

LSS394703	Sling - Liko Original H/B MOD 200 - Plastic Coated Net - w/belt - Sml - Child	\$ 425.00	\$ 212.50	2
LSS394730	Sling - Liko ComfortSling - Mesh - XL 300kg	\$ 599.00	\$ 299.50	4
LSS394750	Sling - Solo Highback Sling - Liko - Med - 2x10 pcs - No Wash - 200kg	\$ 1,629.00	\$ 814.50	1
LSS394760	Sling - Solo Highback Sling - Liko - Lge - 2x10 pcs - No Wash - 200kg	\$ 1,629.00	\$ 814.50	1
MBT401105	Battery - Topin - 6 Volt 5Amp/Hour - for Weigh Chairs	\$ 31.10	\$ 9.33	1
MCU403116	Catheter Bag Holder - Single Unit	\$ 28.90	\$ 8.67	4
MDC411070	Dressing MEPIFORM 4x30cm	\$ 43.70	\$ 13.11	3
MDC414860	Bag ULTIMATE "T" BOX/100	\$ 392.00	\$ 117.60	10
MWS442350	Wheelchair Folding - Aspire METRO - Att Prop - MIDNIGHT BLUE - 100kg ^	\$ 529.00	\$ 264.50	1
MWS442370	Wheelchair Folding - Aspire METRO X - Self Prop - METALLIC ORANGE - 100kg ^	\$ 679.00	\$ 339.50	1
MWS449205	Wheelchair Folding Light - Action 3NG - Self Propelled - 16in (405mm) - 24in solid tyres	\$ 1,329.00	\$ 664.50	1
MWS449900	Wheelchair - Aspire Ward Transfer Chair - 470mm - 140kg	\$ 995.00	\$ 497.50	1
MWS450300	Meal Tray - Aspire Swish Ward Transport Chairs	\$ 495.00	\$ 346.50	4
PED490110	Mattress - Malvestio Smile Paediatric Mattress - 1580x790x140mm - 150kg	\$ 749.00	\$ 524.30	1
PED490550	Seating Support - Firefly - GoTo Size 1 w/Advanced Headrest- Floorsitter- Cushion - Red Vinyl	\$ 1,769.00	\$ 1,238.30	1
PED490615	Seating Support - Firefly - GoTo Size 2 Headrest & Pressure Cushion - Blue	\$ 871.70	\$ 610.19	1
PTA484300	Stockings VENOSAN 4001 Thigh Lace Top AGH C.T X-Large	\$ 170.00	\$ 119.00	1
PTA492600	HipSaver Open Bottom Veterans Kit Medium	\$ 539.00	\$ 377.30	4
PTA493400	HipSaver - SlimFit - High Compliance - Female - XXL	\$ 141.00	\$ 98.70	1
PTA493705	HipSaver - SlimFit - High Compliance - Male - X Sml	\$ 141.00	\$ 98.70	1
PTA494061	HipSaver - QuickChange - Pant Only - Male - XXL	\$ 65.00	\$ 45.50	1
PTA494065	HipSaver - Slim Fit - Pant Only - Male - XL	\$ 65.00	\$ 45.50	1
PTA505100	Hip Protector Pants Male Large (S.M.S.)	\$ 37.00	\$ 25.90	1
PTA505300	Hip Protector Pads (S.M.S.)	\$ 46.00	\$ 32.20	1
PTA507449	Cushion - Positioning - Askle Circular Cushion VCP12	\$ 322.00	\$ 225.40	1
PTA509350	Heel Protector - Maxxcare Heel Pro Evolution - Small	\$ 199.00	\$ 139.30	2
PTA509355	Heel Protector - Maxxcare Heel Pro Evolution - Medium	\$ 199.00	\$ 139.30	12
PTA509360	Heel Protector - Maxxcare Heel Pro Evolution - X Large	\$ 199.00	\$ 139.30	10
PTA511000	HipSaver - Nursing Home - Starter Kit (3-Pants, 1-Pads) - Male - Sml	\$ 269.00	\$ 188.30	1
PTA512060	HipSaver - Nursing Home - Pant Only - Male - XXL	\$ 65.00	\$ 45.50	2
PTA512520	HipSaver - Open Bottom - Pant Only - Lge	\$ 65.00	\$ 45.50	7
PTA512530	HipSaver - Open Bottom - Pant Only - Med	\$ 65.00	\$ 45.50	2
PTA512540	HipSaver - Open Bottom - Pant Only - Sml	\$ 65.00	\$ 45.50	1
PTA512947	HipSaver QC Veterans Kit Female Large	\$ 529.00	\$ 370.30	1
PTA513797	HipSaver - SlimFit - High Compliance w/ Tailbone Pad - Female - X Sml	\$ 170.00	\$ 119.00	1
PTA513798	HipSaver - SlimFit - High Compliance w/ Tailbone Pad - Female - Sml	\$ 170.00	\$ 119.00	1
PTA513800	HipSaver - SlimFit - High Compliance w/ Tailbone Pad - Female - Lge	\$ 170.00	\$ 119.00	1
PTA513801	HipSaver - SlimFit - High Compliance w/ Tailbone Pad - Female - XL	\$ 170.00	\$ 119.00	1
PTA513802	HipSaver - SlimFit - High Compliance w/ Tailbone Pad - Female - XXL	\$ 170.00	\$ 119.00	1
PTA513803	HipSaver - SlimFit - High Compliance w/ Tailbone Pad - Male - X Sml	\$ 170.00	\$ 119.00	1
PTA513804	HipSaver - SlimFit - High Compliance w/ Tailbone Pad - Male - Sml	\$ 170.00	\$ 119.00	1
PTA513807	HipSaver - SlimFit - High Compliance w/ Tailbone Pad - Male - XL	\$ 170.00	\$ 119.00	1
PTA513808	HipSaver - SlimFit - High Compliance w/ Tailbone Pad - Male - XXL	\$ 170.00	\$ 119.00	1
PTA514410	Seat Adjustable - Jay - 15-18in - Solid with Lip Kit	\$ 559.00	\$ 391.30	1
PTA514800	Universal Post Bracket - Jay - Bottom	\$ 156.00	\$ 109.20	3
PTA514900	Universal Post Bracket - Jay - Top - Quick Release	\$ 87.00	\$ 60.90	3
PTA515400	Build Up Abductor - Jay - 3in - for Jay Cushions - GS	\$ 116.00	\$ 81.20	1
PTA520499	Boot - Sheepskin Hi Temp - Tall Slipper - Open Toe Small	\$ 276.00	\$ 193.20	1
PTA522700	Mattress Pressure Low Level & Pump	\$ 359.00	\$ 251.30	1
PTA524800	Stockings JOBST Ultrasheer Knee High Closed Toe 15-20mm Small Sun Brown(S.M.S.)	\$ 52.00	\$ 36.40	6
PTA532964	DermaSaver Shin Tube Small	\$ 77.00	\$ 53.90	1
PTA533010	Limb Protection - DermaSaver - Arm Tube with Double Elbow - X Sml	\$ 93.00	\$ 65.10	1
PTA533324	Limb Protection - GelBodies - Heel & Elbow Protector - Sml	\$ 118.00	\$ 82.60	2
PTA534101	Limb/Cast Protection - Limbo - Adult Hand - Waterproof	\$ 51.00	\$ 35.70	1
PTA534102L	Limb/Cast Protection - Limbo - Adult Below Elbow - Waterproof - Lge	\$ 61.00	\$ 42.70	1
PTA534111S	HipSaver - Track Pant w/ Knee Pad - Regular - Sml	\$ 181.00	\$ 126.70	1
PTA534112M	HipSaver - Track Pant w/ Tailbone & Knee Pad - Regular - Med	\$ 208.00	\$ 145.60	1
PTA534230	Body Wedge - Snoooooze - Wedge PAEDIATRIC	\$ 210.43	\$ 147.30	1
PTB532180	WChair Back - JAY J2 - Deep Contour 14in	\$ 1,429.00	\$ 714.50	1
PTB532310	WChair Back - JAY J2 - 16in	\$ 989.00	\$ 494.50	1
PTB532330	WChair Back - JAY J2 - 17in	\$ 989.00	\$ 494.50	2
PTB532340	WChair Back - JAY J2 - 17in Tall	\$ 1,129.00	\$ 564.50	1
PTB532444	Backrest Icon 16" Tall	\$ 1,439.00	\$ 719.50	1
PTB533685	J3 Back Posterior & Deep Contour Shoulder Height Small 18in	\$ 1,529.00	\$ 764.50	1
PTB533823	WChair Back - JAY J3 - Posterior/Lateral Pelvic	\$ 86.00	\$ 43.00	2
PTB533837	WChair Back - JAY J3 - Laterals - Swing Away 6x5in	\$ 759.00	\$ 379.50	1
PTB534500	Matrix Backrest Cover Size 4	\$ 589.00	\$ 294.50	1
PTB535174	Comfi Back Support - 18in Wide x 18in High with tear-off comfort cubes	\$ 1,619.00	\$ 809.50	1
PTB535292	Spex Adjustable Back Support 15in Wide x 18in High	\$ 2,089.00	\$ 1,044.50	1
PTB538850	WChair Back - Dreamline Contour Backrest - Fixed - 18in x 22in	\$ 1,449.00	\$ 724.50	1
PTB541420	Wheelchair Backrest - WheelAir - SlingBack - Small - 360mm	\$ 2,359.00	\$ 707.70	1
PTB541480	Wheelchair Backrest - WheelAir - SlingBack - Medium - 360mm	\$ 2,489.00	\$ 746.70	1
PTC536160	J3 Harness Centre Opening with Buckle Large	\$ 398.00	\$ 119.40	1
PTC542150	Cushion - Seat Pad Silicone Fibre 16x16	\$ 44.60	\$ 13.38	1
PTC543500	Cushion Diffuser Water Proof Overslip Cover (S.M.S.)	\$ 31.00	\$ 9.30	1
PTC546400	Cushion - ROHO Single Compartment HP 16.75x9.25in (9x5cell)	\$ 1,000.00	\$ 500.00	1

WAREHOUSE CLEARANCE

ALL PRODUCTS BETWEEN **50-70% OFF**

HURRY, WHILST STOCK LASTS!

PTC548700	Cushion - Spenco WChair -Full Padding & Tiedown	\$ 235.00	\$ 117.50	1
PTC549900	Cushion - Cloud - Otto Bock - 19x18in	\$ 1,279.00	\$ 639.50	1
PTC551580	Cushion Vicair Adjuster 10 & incotec Cover - 16x16in	\$ 899.00	\$ 449.50	2
PTC554200	Cushion Infinity Flow Gel 18 x 16" (S.M.S.)	\$ 754.00	\$ 377.00	1
PTC556300	Cushion Cover - ROHO Waterproof LP 15x15in	\$ 179.00	\$ 53.70	1
PTC556600	Cushion Cover - ROHO Waterproof HP 16.75x18.5in	\$ 179.00	\$ 53.70	2
PTC559000	Cushion - JAY - GS 12x15in	\$ 660.00	\$ 330.00	1
PTC561271	Cushion - Jay X2 16 x18in	\$ 759.00	\$ 379.50	2
PTC561274	Cushion - Jay X2 - 17 x17in	\$ 759.00	\$ 379.50	1
PTC563100	Cushion - JAY J2 - Recline 15.5x20in	\$ 1,159.00	\$ 579.50	2
PTC569300	Cushion - JAY - Duo 16x18in	\$ 989.00	\$ 494.50	1
PTC570000	JAY Cushion Xtreme 15x16" Tritex (S.M.S.)	\$ 647.00	\$ 323.50	1
PTC571800	Cushion - JAY - Combi 16x16in	\$ 256.00	\$ 128.00	1
PTC571980	Cushion Cover - JAY - Lite 16x16in	\$ 145.00	\$ 43.50	7
PTC573800	Cushion Cover - JAY - Stretch 15.5x18in	\$ 87.00	\$ 26.10	1
PTC576200	Cushion Cover - JAY J2 - 18x20in - Recline Stretch	\$ 156.00	\$ 46.80	1
PTC576500	Cushion Cover - JAY J2 - 15.5x16in - Recline Air Exch	\$ 126.00	\$ 37.80	1
PTC577500	Cushion Cover - JAY J2 - 17x17in - Recline Stretch	\$ 126.00	\$ 37.80	1
PTC582100	JAY Fluid Pad Deep Contour 20x20in	\$ 1,099.00	\$ 549.50	1
PTC590400	Cushion - ROHO Quadtro Select HP 12x13in (6x7cell)	\$ 1,119.00	\$ 559.50	1
PTC594920	Cushion - ROHO Hybrid Elite 1R 16x18in Sensor Ready	\$ 1,079.00	\$ 539.50	1
PTC594920SR	Cushion - ROHO Hybrid Elite 1R 16x18in w/HD cover Sensor Ready	\$ 1,079.00	\$ 539.50	3
PTC594921	Cushion - ROHO Hybrid Elite 1R 16x20in Sensor Ready	\$ 1,079.00	\$ 539.50	1
PTC594929	Cushion - ROHO Hybrid Elite 1R 17x16in Sensor Ready	\$ 1,079.00	\$ 539.50	1
PTC594935	Cushion - ROHO Hybrid Elite 1R 18x16in Sensor Ready	\$ 1,079.00	\$ 539.50	1
PTC594945	Cushion - ROHO Hybrid Elite 1R 18x20in	\$ 839.00	\$ 419.50	1
PTC594956	Cushion - ROHO Hybrid Elite 2R 15x15in	\$ 1,079.00	\$ 539.50	1
PTC596150	Cushion Cover - ROHO Universal HP 18.5x16in	\$ 153.00	\$ 45.90	1
PTC596500	Cushion Cover - ROHO Standard HP 18.5x22in	\$ 153.00	\$ 45.90	1
PTC600500	Cushion Synergy Structure 18 x 16"	\$ 662.00	\$ 331.00	1
PTC600600	Cushion Synergy Structure 18 x 18"	\$ 772.00	\$ 386.00	1
PTC601450	Cushion - Synergy Spectrum - Air Contour - w/ R/H valve 18 x 18in	\$ 836.20	\$ 418.10	1
PTC604100	Cushion Cover Keyhole Cushion Overslip (S.M.S.)	\$ 16.00	\$ 8.00	1
PTC609000	Cushion - Vicair - O2 Adjuster 6 Top cover - 18x16inch	\$ 819.00	\$ 409.50	2
PTC609120	Cushion - Vicair - O2 Adjuster 6 Incotec - 18x16inch	\$ 819.00	\$ 409.50	1
PTC609170	Cushion - Vicair - O2 Adjuster 6 Incotec - 20x20inch	\$ 819.00	\$ 409.50	4
PTC609370	Cushion - Vicair - Liberty Profile Air Exchange - 18x18inch	\$ 350.00	\$ 175.00	3
PTC609490	Cushion - Vicair - O2 Active 9 - 20x16inch	\$ 849.00	\$ 424.50	2
PTC609500	Cushion - Vicair - O2 Active 10 - 20x18inch	\$ 849.00	\$ 424.50	1
PTC609800	Cushion - Vicair - Adjuster 12 O2 Incotec - 16x18inch	\$ 999.00	\$ 499.50	1
PTC620630	Cushion - Dreamline Support Cushion - 18in x 18in	\$ 929.00	\$ 464.50	1
PTC620900	Cushion - Dreamline Contour Cushion - 16in x 18in	\$ 979.00	\$ 489.50	1
PTC621975	Cushion - Axiom - Standard - 18 x 16	\$ 157.00	\$ 78.50	1
PTC622510	Cushion - Pride Stealth Spectrum Cushion - Foam - 16in x 18in	\$ 467.00	\$ 233.50	1
PTM636560	Mattress Overlay - Ultimate 75 - Alt Air - 21x8cm cell - 2000x880x80mm	\$ 1,184.50	\$ 592.25	2
PTM636565	Mattress Cover - Concave Positioning Cover - Alt Air - 2000x880x100	\$ 889.00	\$ 266.70	1
PTM636625	AirLo - Thermal Hyperstretch Microfibre Coverlet	\$ 494.00	\$ 148.20	3
PTM638015	Mattress - Aspire HYBRID - King Single - 9 Cell - 1980x1050x170mm - 180kg (Pump Sep: PTM638030)	\$ 1,939.00	\$ 969.50	1
PTM638020	Mattress - Aspire HYBRID ACUTE - Single - 9 Cell - 1980x895x170mm - 230kg (Pump Sep: PTM638030)	\$ 1,939.00	\$ 969.50	2
PTM639030	Mattress Cover - Aspire Concave Positioning Cover - 2000x900x125mm Active Air 5 Single	\$ 519.00	\$ 155.70	1
PTM639300T	Mattress Replacement - Aspire ACTIVE AIR 8 TURN - Single - 20 Cells - 2000x900x200mm - 200kg	\$ 2,821.00	\$ 1,410.50	1
PTM639400	Mattress Replacement - Aspire ACTIVE AIR ACUTE 8 - Single - 24 Cells - 2000x900x200mm - 250kg	\$ 3,775.00	\$ 1,887.50	1
PTS640021	JOBSTTravel Socks Black Size 1	\$ 58.00	\$ 17.40	1
PTS640022	JOBSTTravel Socks Black Size 2	\$ 58.00	\$ 17.40	2
PTS640024	JOBSTTravel Socks Black Size 4	\$ 58.00	\$ 17.40	2
PTS640040	Stockings - JOBST for Men - Knee - 15-20mmHg - C/Toe - Sml - Blk	\$ 70.00	\$ 21.00	4
PTS640080	Stockings - JOBST for Men - Knee - 15-20mmHg - C/Toe - Sml - Navy 1pair	\$ 62.00	\$ 18.60	1
PTS640270	Stockings - JOBST UltraSheer - Knee - 15-20mmHg - C/Toe - Lge - Suntan	\$ 49.00	\$ 14.70	1
PTS640315	Stockings - JOBST UltraSheer - Thigh w/ Grip Top - 15-20mmHg - C/Toe - XL - Blk	\$ 97.00	\$ 29.10	2
PTS640420	JOBST UltraShr Waist 15-20 Navy Small (S.M.S.)	\$ 86.00	\$ 25.80	1
PTS640500	JOBST UltraShr Mater 15-20 Natural Medium	\$ 155.00	\$ 46.50	1
PTS640580	Stockings - JOBST for Men - Knee - 20-30mmHg - C/Toe - Lge - Navy	\$ 93.00	\$ 27.90	1
PTS640640	Stockings - JOBST UltraSheer - Knee - 20-30mmHg - C/Toe - XL - Blk	\$ 103.00	\$ 30.90	1
PTS640690	Stockings - JOBST UltraSheer - Knee - 20-30mmHg - C/Toe - Med - Suntan	\$ 93.00	\$ 27.90	3
PTS640740	Stockings - JOBST UltraSheer - Thigh w/ Grip Top - 20-30mmHg - C/Toe - Sml - Natural	\$ 128.00	\$ 38.40	2
PTS640781	Stockings - JOBST UltraSheer - Waist - 20-30mmHg - C/Toe - Lge - Blk	\$ 207.00	\$ 62.10	1
PTS640800	Stockings - JOBST UltraSheer - Waist - 20-30mmHg - C/Toe - Lge - Natural	\$ 202.00	\$ 60.60	1
PTS640815	Stockings - JOBST UltraSheer - Waist - 20-30mmHg - C/Toe - Lge - Suntan	\$ 181.00	\$ 54.30	1
PTS640830	Stockings - JOBST for Men - Knee - 30-40mmHg - C/Toe - Sml - Blk	\$ 113.00	\$ 33.90	2
PTS640920	Stockings - JOBST UltraSheer - Knee - 30-40mmHg - C/Toe - Med - Natural	\$ 102.00	\$ 30.60	2
PTS640940	Stockings - JOBST UltraSheer - Knee - 30-40mmHg - C/Toe - XL - Natural	\$ 102.00	\$ 30.60	1
PTS641210	Stockings - JOBST Relief - Knee - 30-40mmHg - Unisex - C/Toe - Lge	\$ 102.00	\$ 30.60	1
PTS641265	Stockings - Jobst Relief - Thigh w/Grip Top - 20-00mmHg - Unisex - C/Toe - Sml	\$ 130.00	\$ 39.00	1
PTS641335	Stockings - JOBST Relief - Waist - 20-30mmHg - P/hose - Unisex - C/Toe - XL	\$ 181.00	\$ 54.30	3
PTS641365	Stockings - JOBST Relief - Thigh w/ Belt (Chaps) Right - 20-30mmHg - Unisex - O/Toe - Lge	\$ 201.00	\$ 60.30	1
PTS641520	TED Stocking Knee Closed Toe Large Long Beige	\$ 53.00	\$ 15.90	2

WAREHOUSE CLEARANCE

ALL PRODUCTS BETWEEN **50-70% OFF**

HURRY, WHILST STOCK LASTS!

PTS641550	Compress Stockings - TED - Thigh Closed Toe - Medium Long	\$ 58.00	\$ 17.40	1
PTS641632	Stockings - TED Continuing Care - Knee - Anti-embolic - C/Toe - White - Lge	\$ 63.00	\$ 18.90	2
PTS641740	Compress Stockings - TED - Knee High Small Long Beige	\$ 42.70	\$ 12.81	3
PTS641820	TED Stocking Thigh High Extra Large & Belt Long (S.M.S.)	\$ 42.00	\$ 12.60	1
PTS641860	Stockings - GloriaMed - Unisex Knee - 20-30mmHg - O/Toe - X Sml - Bge	\$ 70.00	\$ 21.00	2
PTS641957	Stockings - Venosan 4002 (Nylon) - Knee - 23-32mmHg - Unisex - O/Toe - Sml - Marokko	\$ 143.00	\$ 42.90	1
PTS641960PTMAS	Stockings - Venosan 4002 (Nylon) - Thigh / Plain - 23-32mmHg - Unisex - C/Toe - Sml - Marokko	\$ 224.00	\$ 67.20	1
PTS641974	Stockings - Venosan 6001 (Nyl/Lyc) - Knee - 18-22mmHg - Unisex - O/Toe - Short - Sml - Blk	\$ 107.00	\$ 32.10	2
PTS641974OLBES	Stockings - Venosan 6001 (Nyl/Lyc) - Knee - 18-22mmHg - Unisex - O/Toe - Long - Sml - Bge	\$ 107.00	\$ 32.10	4
PTS641974OSBES	Stockings - Venosan 6001 (Nyl/Lyc) - Knee - 18-22mmHg - Unisex - O/Toe - Short - Sml - Bge	\$ 107.00	\$ 32.10	1
PTS641974OSBLXL	Stockings - Venosan 6001 (Nyl/Lyc) - Knee - 18-22mmHg - Unisex - O/Toe - Short - XL - Blk	\$ 107.00	\$ 32.10	1
PTS642001	Stockings - Mediven Elegance - Knee - 20-30mmHg - C/Toe - Short - (Sz I) - Bge	\$ 130.00	\$ 39.00	1
PTS642017	Stockings - Mediven Comfort - Knee (Short) - 20-30mmHg - C/Toe - (Sz I) - Bge	\$ 136.00	\$ 40.80	2
PTS648485	Stockings - Sigvaris Cotton - Thigh w/Grip Top - 18-21mmHg L Norm Nature Closed Toe	\$ 278.00	\$ 83.40	1
PTS650390	Stockings - Venosan 4001 (Nylon) - Thigh w/ Plain Top - 18-22mmHg - Unisex - O/Toe - XLge - Long - M	\$ 219.00	\$ 65.70	1
PTS650410	Stockings - Venosan 4001 (Nylon) - Waist - 18-22mmHg - Unisex - P/hose - C/Toe - XLge - Long - Mexic	\$ 255.00	\$ 76.50	1
PWS655320	Power Wheelchair - Whill - Model Ci - (Black Armrests) - 450mm - 114kg	\$ 7,489.00	\$ 3,744.50	1
SCA668220	Scooter Acc - Cane Holder Double (Rear) w/flag holder (requires L Piece)	\$ 115.00	\$ 34.50	4
SCT670540	Scooter - Pride - Celebrity Deluxe XL - 4 Wheel Red (with Rear bag + Cane Holder)	\$ 6,379.00	\$ 3,189.50	1
TSE672620	Gloves Multiglide PER PAIR (S.M.S.)	\$ 37.00	\$ 11.10	2
TSE672870	Transfer - Walk Belt - Hemco - Large	\$ 45.92	\$ 13.78	17
TSE672890	Transfer - Walk Belt - Hemco - Small	\$ 45.92	\$ 13.78	17
TSE674400	Slide Sheet - SlipperySally - Single Patient Use - 2.0m x 1.45m (Box of 50)	\$ 699.00	\$ 209.70	2
TSE674410	Slide Sheet - SlipperySally - Single Patient Use (BAGGED) - 2.0m x 1.45m (Box of 50)	\$ 719.00	\$ 215.70	1
TSE674420	Repositioning Sheet - Sally Sheet - Single Patient Use - Bariatric 1200mm wide - 454kg (Box of 5)	\$ 474.00	\$ 142.20	1
TSE674430	Repositioning Sheet - Sally Sheet - Single Patient Use - Bariatric 1500mm wide - 454kg (Box of 5)	\$ 569.00	\$ 170.70	1
TSE674440	Air Assisted Lateral Transfer Mat Single Patient Use - 860 W x 1980 L (mm) - (Box of 10)	\$ 2,499.00	\$ 749.70	3
TSE674450	Air Assisted Lateral Transfer Mat - Air Supply 1200W	\$ 2,459.00	\$ 737.70	1
TSE674460	Air Assisted Lateral Transfer Mat - 1.7m Replacement Supply Hose	\$ 445.00	\$ 133.50	1
TSE674470	Air Assisted Lateral Transfer Mat - Transport Stand - Single Basket	\$ 1,439.00	\$ 431.70	1
TSE674480	Air Assisted Lateral Transfer Mat - Disposable Waterproof Flat Sheet (Box of 50)	\$ 300.00	\$ 90.00	2
UPA200210	DDC Detergent Pump - MAXI, Optima 3	\$ 1,239.00	\$ 371.70	1
UPA200420	Temperature Monitoring Kit - Maxi & Optima 3	\$ 1,679.00	\$ 503.70	1
WAA686225	Crutch Padding Grey Hand Grip - Red Dot	\$ 10.70	\$ 3.21	62
WAA686230	Crutch Padding Grey Hand Grip - Bariatric Crutches - Wing Nut	\$ 9.46	\$ 2.84	120
WAA686250	Crutch Acc - Aspire - Grey Padded Hand Grip	\$ 20.70	\$ 6.21	63
WAA686370	Tip Ferrule Bowlers 16mm for Walking Sticks	\$ 36.80	\$ 11.04	1
WAA686375	Tip Ferrule Bowlers 19mm for Walking Sticks	\$ 36.80	\$ 11.04	2
WAA686911	Aspire Walker Accessories - Underseat Basket - XL	\$ 31.80	\$ 9.54	8
WAA686925	Walking Stick Holder - Aspire Walkers - Superlite/ XL	\$ 14.60	\$ 4.38	2
WAA686980	Walker Acc - Gutter Walker - Forearm Supports Cuff Pad	\$ 118.00	\$ 35.40	2
WAA695615	Crutch & Cane Holder - Aspire VOGUE INDOOR Walkers	\$ 12.50	\$ 3.75	40
WAC692495	Crutches Underarm - Alum - Wing Nut - 110kg - Adult Tall	\$ 38.00	\$ 11.40	33
WAF699216	Accessory Pack (IV Pole, Oxy Holder & Storage Bag) - Freedom Titan Forearm Walker	\$ 92.00	\$ 27.60	3
WAF705600BL	Seat Walker - Aspire PUSH DOWN - Blue - 7in Wheels - 130kg	\$ 207.00	\$ 103.50	41
WAF705600RE	Seat Walker - Aspire PUSH DOWN - Red - 7in Wheels - 130kg	\$ 207.00	\$ 103.50	15
WAF705900	Seat Walker - Envoy - H-Brakes - Heavy Duty 170kg - 8in Wheels	\$ 229.00	\$ 68.70	1
WAF709300	Forearm Walker - K Care - Walking Tutor - Junior w/ Glides - 670-960mm	\$ 1,509.00	\$ 452.70	1
WAF709650CH	Seat Walker - Aspire SUPERLITE Walker - Champagne - HBrakes - 6in Wheels - 113kg ^	\$ 219.00	\$ 109.50	1
WAF750150	Seat Walker - Freedom - Stride - Red - 135kg	\$ 369.00	\$ 184.50	105
WAF751100	Seat Walker - Saljol - Page Indoor Walker - Small - Charcoal - 150kg	\$ 380.00	\$ 114.00	3
WAF751110	Seat Walker - Saljol - Page Indoor Walker - Medium - Charcoal - 150kg	\$ 380.00	\$ 114.00	3
WAS744710	Walking Stick - Carequip - Ht Adj - T Handle - Folding - Alum	\$ 73.00	\$ 21.90	1
WAS748110	Walking Stick - Freedom - Cobra Cane - 710-940mm - 110kg	\$ 56.00	\$ 16.80	3
WAS748140	Walking Stick - Freedom - Crook Handle - Black - 770-970mm - 110kg	\$ 28.50	\$ 8.55	1
WAS748160	Walking Stick - Freedom - Crook Handle - Bronze - 770-970mm - 110kg	\$ 28.50	\$ 8.55	2
WAS748180	Walking Stick - Freedom - Crook Handle - Floral - 770-970mm - 110kg	\$ 28.50	\$ 8.55	5
WAS748190	Walking Stick - Freedom - Crook Handle - Leaves - 770-970mm - 110kg	\$ 28.50	\$ 8.55	12
WAS748210	Walking Stick - Freedom - Crook Handle - Tartan - 770-970mm - 110kg	\$ 28.50	\$ 8.55	3
WAS748230	Walking Stick - Freedom - Folding - Black - 800-900mm - 114kg	\$ 28.50	\$ 8.55	1
WAS748280	Walking Stick - Freedom - Offset Handle - Tartan - 770-970mm - 110kg	\$ 24.00	\$ 7.20	2
WAS748340	Rubber Tip for Quad Cane - 12mm - Freedom	\$ 5.40	\$ 1.62	6
WAS748370	Walking Stick - Freedom - Swan Neck - Bronze - 770-970mm - 110kg	\$ 28.50	\$ 8.55	1
WAS748380	Walking Stick - Freedom - Swan Neck - Chrome - 770-970mm - 110kg	\$ 28.50	\$ 8.55	3
WAS748390	Walking Stick - Freedom - Crook Handle - Chrome - 770-970mm - 110kg	\$ 28.50	\$ 8.55	2
WAS748400	Walking Stick - Freedom - Swan Neck - Floral - 770-970mm - 110kg	\$ 28.50	\$ 8.55	4
WAS748410	Walking Stick - Freedom - Swan Neck - Leaves - 770-970mm - 110kg	\$ 28.50	\$ 8.55	1
WAS748420	Walking Stick - Freedom - Swan Neck - Tartan - 770-970mm - 110kg	\$ 28.50	\$ 8.55	1
WAS748510	Walking Stick - Freedom - T Handle Ergo Plastic - Floral - 770-970mm - 110kg	\$ 28.50	\$ 8.55	2
WAS748540	Walking Stick - Freedom - T Handle Plastic - Bronze - 770-970mm - 110kg	\$ 28.50	\$ 8.55	1
WAS748580	Walking Stick - Freedom - T Handle Plastic - Tartan - 770-970mm - 110kg	\$ 28.50	\$ 8.55	2
WAT758800	Bedside Cabinet - Mini Duke - 3 Drawer w/Lock - Sublime Teak Laminate	\$ 633.08	\$ 189.92	3
WAT760800	Wash Bowl - Green Autoclavable - 240mm x 110mm - 2000ml	\$ 21.40	\$ 6.42	8
WAT761400B	Bedside Cabinet - 2 Drawers + 1 Dump Drawer With Lock (Beech)	\$ 378.50	\$ 113.55	1
WAT774600	IV Stand Deluxe for wheelchairs	\$ 182.00	\$ 54.60	2
WAT805400	Instrument Trolley - Kerry - 500x500 - 1 Drawer S/S - Epoxy PC Frame	\$ 580.00	\$ 290.00	3

WAREHOUSE CLEARANCE

ALL PRODUCTS BETWEEN **50-70% OFF**

HURRY, WHILST STOCK LASTS!

WAT805700	Instrument Trolley - Kerry - 500x500 - 1 Drawer S/S - Chrome Frame	\$ 1,019.00	\$ 509.50	1
WAT807800	Dressing Trolley - Kerry - 500x500 - Flat Top - Stainless Steel	\$ 985.60	\$ 492.80	2
WAT809100	Medicine Cart - Kerry - Webster Cart - 3 Drawers (48) - Beige	\$ 3,223.80	\$ 1,611.90	1
WAT816650	FlexiCart Options - Folding Shelf Heavy Duty with Flat Top	\$ 719.00	\$ 215.70	1
WAT820180	IV Pole - Aspire - 5 Star Polyurethane Base	\$ 188.00	\$ 94.00	5
WAT828200	Bedside Cabinet - Kerry - 1 Drawer w/ Cupboard - Economy (Custom Code - Price on Request)	\$ 739.00	\$ 369.50	2
WAT900635	Dressing Trolley - Aidacare - S/Steel - 490x490x900 - 2 Drawer	\$ 1,059.00	\$ 529.50	1
WAT900735	Emergency Cart - Aidacare - 650x480x980mm - Red - 5 Drawers	\$ 3,229.00	\$ 1,614.50	2
WAT905110	Emerg Cart TRIBUTE 90 PKG - 8H DR311 - BRK AWY Lock & RHS Handle - CPR Bd, IV, CO2, Defib, Waste - R	\$ 3,579.00	\$ 1,789.50	1
WAT905350-O	Isolation Cart Avalo 9H - DR003 w/ No Lock & RHS Handle - YELLOW (dividers incl)	\$ 2,374.00	\$ 1,187.00	1
WAT907185	Avalo Multi-Cavity Stotage - Complete	\$ 237.00	\$ 71.10	1
WAT907190	Cart Accessory - Aidacare/Tribute - Anesthesia Storage Tower Assy (needs 2 x A&E Mount Brackets 1214	\$ 180.00	\$ 54.00	1
WAT907250	Avalo Handle & End Panel- 10H (K) - Light Cream/Grey	\$ 425.00	\$ 127.50	4
WAT907252	Avalo Handle & End Panel- 10H (K) - Ext Red	\$ 425.00	\$ 127.50	1
WAT907254	Avalo Handle & End Panel- 10H (K) - Light Blue	\$ 425.00	\$ 127.50	1
WAT907260	Avalo Handle & End Panel- 9H (K) - Light Cream/Grey	\$ 425.00	\$ 127.50	2
WAT907264	Avalo Handle & End Panel- 9H (K) - Light Blue	\$ 425.00	\$ 127.50	2
WAT907280	Avalo End Panel only - 10H (K) - Light Cream/Grey	\$ 263.00	\$ 78.90	3
WAT907325	Cart Accessory - Aidacare/Tribute - Draw Tab Filler - Grey	\$ 2.31	\$ 0.69	1
WAT907335	Drawer Divider - Tall - Grey (6 & 9in Drawers) - TRIBUTE	\$ 26.60	\$ 7.98	1
WAT907460	Avalo Drawer Slide Assembly (K)	\$ 107.00	\$ 32.10	16
WAT907470	Avalo Mechanical Lock Assembly (K)	\$ 336.00	\$ 100.80	1
WAT907520	Avalo Keyless Pcb	\$ 390.00	\$ 117.00	1
WAT907570	Avalo Lock Tab Kit (K)	\$ 71.00	\$ 21.30	1
WAT907600	Avalo Lock Bezel Override Cover	\$ 13.80	\$ 4.14	1
WAT907620	Avalo Lock Bezel Exterior Filler - Light Cream/Grey	\$ 43.90	\$ 13.17	1
WAT907665	Avalo Writing Surface - Light Cream	\$ 173.00	\$ 51.90	4
WAT907682	Avalo Bumper - Specify Colour - Ext Blue	\$ 98.00	\$ 29.40	2
WAT907716	Avalo Lock Tab 3.4 Repeat	\$ 11.90	\$ 3.57	3
WAT907717	Avalo Lock Tab 3.4 Poly	\$ 13.10	\$ 3.93	4
WCA848170	WChair Elevating Legrest - Advanced - Right - 320-410mm - Aspire Assist2/Transit2/Evoke2/Rehab RS&R	\$ 183.00	\$ 54.90	1
WCA848190	WChair Elevating Legrest - Advanced - Right - 450-510mm - Aspire Assist2/Transit2/Evoke2/Rehab RS&R	\$ 183.00	\$ 54.90	1
WCA848210	WChair Elevating Legrest - Articulating - Right - 570/610mm - Aspire Assist2/Evoke2	\$ 183.00	\$ 54.90	1
WCA848325	WChair Contoured Padded Overlay - 300mm - Aspire Evoke2 JNR	\$ 99.00	\$ 29.70	5
WCA848330	WChair Contoured Padded Overlay - 350mm - Aspire Evoke2	\$ 99.00	\$ 29.70	6
WCA848490	WChair Camber Kit - Aspire Evoke2	\$ 17.40	\$ 5.22	2
WCA848510	WChair Drum Brakes Kit - 16in Wheels - Aspire Rehab RS/Rehab RX	\$ 326.00	\$ 97.80	1
WCA848520	WChair Drum Brakes Kit - 20in Wheels - Aspire Rehab RS/Rehab RX	\$ 341.00	\$ 102.30	3
WCA848525	WChair Drum Brakes Kit - 24in Wheels - Aspire Rehab RS/Rehab RX/Rehab RX JNR	\$ 349.00	\$ 104.70	11
WCA848530	WChair Depth Extension Kit - 50mm - Black - Aspire Rehab RS/Rehab RX	\$ 25.80	\$ 7.74	6
WCA848540	WChair Depth Extension Kit - 100mm - Black - Aspire Rehab RS/Rehab RX	\$ 32.30	\$ 9.69	10
WCA848550	WChair Therapist Kit for Hand Adjustment - Aspire Rehab RS/Rehab RX	\$ 38.30	\$ 11.49	5
WCA848590	WChair Extension/Stroller Handle - Aspire Rehab RX	\$ 311.00	\$ 93.30	2
WCA858005	WChair Elevating L/rests- alu f/plates - Heavy Duty WChair (20-24in)	\$ 152.00	\$ 45.60	2
WCA858203	90" Lap Belt	\$ 49.90	\$ 14.97	2
WCA858206	WChair Belt - 70in Lap Belt	\$ 58.00	\$ 17.40	2
WCA858390	WChair Elevating Legrest - Aspire EVOKE 400/450/500/HD - Black Right	\$ 120.00	\$ 36.00	1
WCA858400	WChair Elevating Legrest - Aspire ASSIST/TRANSIT - Black Left	\$ 110.00	\$ 33.00	4
WCA858410	WChair Elevating Legrest - Aspire ASSIST/TRANSIT - Black Right	\$ 110.00	\$ 33.00	2
WCA858450	WChair Stump Support - Aspire EVOKE - Black Right	\$ 105.00	\$ 31.50	4
WCA858460	WChair Anti Tip Bar - Aspire EVOKE	\$ 52.00	\$ 15.60	1
WCA858462	WChair Anti Tip Bar - Aspire ASSIST/TRANSIT	\$ 25.00	\$ 7.50	28
WCA858472	Adjustable Backrest Cane (Pair) - Aspire EVOKE - 19-21in Height	\$ 60.00	\$ 18.00	25
WCA858590	WChair Belt w/ Buckle - Attachment	\$ 93.00	\$ 27.90	5
WCA858610	WChair Belt w/ Plastic Clip Continuous	\$ 93.00	\$ 27.90	5
WCA858760	Combination Head/Neckrest - Ottobock	\$ 246.00	\$ 73.80	2
WCA858823	Headrest Clamp Heavy Duty (New Style) - Black	\$ 268.00	\$ 80.40	2
WCA858870	Side Thoracic Support Mounting Kit Dropped Left Black	\$ 195.00	\$ 58.50	1
WCA858880	Side Thoracic Support Mounting Kit Dropped Right Black	\$ 195.00	\$ 58.50	1
WCA858940	Side Thoracic Support Large	\$ 219.00	\$ 65.70	1
WCA859000	Channel Forearm Pad Small - Ottobock	\$ 239.00	\$ 71.70	6
WCA859015	Side Panel Armrest Pad	\$ 384.00	\$ 115.20	1
WCA859160	Modular Hand Pad Horn Large Left	\$ 125.00	\$ 37.50	2
WCA859250	Hand Support Flat - Ottobock - Size Medium	\$ 125.00	\$ 37.50	1
WCA859340	Clinical Mounting Kit 52-62cm Black	\$ 356.00	\$ 106.80	1
WCA859447	Start Wheel Lock Lever Extension	\$ 33.40	\$ 10.02	1
WCA859770	WChair Oxygen Bottle Holder - Aspire ASSIST/TRANSIT	\$ 186.00	\$ 55.80	4
WCA860241	WChair Pelvic Belt - Push Button - Centre Pull 4 Point - Lge	\$ 301.00	\$ 90.30	2
WCA870100	Power Wheelchair Accessory - Whill C1 - Arm Covers - Black (Left)	\$ 195.00	\$ 58.50	2
WCA870110	Power Wheelchair Accessory - Whill C1 - Arm Covers - Black (Right)	\$ 195.00	\$ 58.50	2
WCA870160	Power Wheelchair Accessory - Whill C1 - Arm Covers - Light Blue (Left)	\$ 195.00	\$ 58.50	1
WCA870240	Power Wheelchair Accessory - Whill C1 - Battery Charger	\$ 225.00	\$ 67.50	2
WCA870300	Power Wheelchair Accessory - Whill C1 - Storage Basket	\$ 246.00	\$ 73.80	1
WCA870320	Power Wheelchair Accessory - Whill C1 - Lithium Ion Battery	\$ 1,549.00	\$ 464.70	4