

# WAREHOUSE CLEARANCE

ALL PRODUCTS BETWEEN **50-70% OFF**

HURRY, WHILST STOCK LASTS!

| Item Number | Item Description  | WAS         | NOW         | QTY |
|-------------|---|-------------|-------------|-----|
| BEA007488   | Alarm - Cura - Cordless Floor Mat - (1220 x 610) includes CU-2619 Transmitter                   | \$ 649.00   | \$ 194.70   | 1   |
| BEA007850   | Alarm - Sofihub TEQ Home (1 x SOFIHUB TEQ Home + 8 Sensors)                                     | \$ 1,419.00 | \$ 709.50   | 1   |
| BEA014605   | Bed Cradle - Aspire - Fixed - Treated Steel ^   | \$ 69.00    | \$ 20.70    | 6   |
| BEA024100   | Bed End Set - CS7 - Amherst 36in w/ Ctrl Panel - Solar Oak                                      | \$ 579.00   | \$ 173.70   | 2   |
| BEA024128   | Assist Rail - Invacare - for Invacare CS7 Wide Deck Bed   | \$ 402.00   | \$ 120.60   | 1   |
| BEA029845   | Bed End Shrouds - Beech - Single - Aspire COMMUNITY CARE Bed (Premium)                          | \$ 300.00   | \$ 90.00    | 1   |
| BEA030050   | Bed Mover/Transporter - Aidacare FL250 Bed (Set of 2))  | \$ 809.00   | \$ 404.50   | 1   |
| BEB046605   | Reprogramming Handset - Safety Egress Height Adjustment - Aidacare FL250/FLX Bed                | \$ 200.00   | \$ 60.00    | 1   |
| BEM047310   | Convuluted Foam Mattress Overlay Cover Protect A Plus PVC 190cmx91cm                            | \$ 38.00    | \$ 19.00    | 1   |
| BTC062500   | Shower Commode Seat - K Care - H/Made - Closed 550mm  | \$ 839.00   | \$ 251.70   | 1   |
| BTC062600   | Shower Commode Seat - K Care - H/Made - Closed 600mm  | \$ 999.00   | \$ 299.70   | 1   |
| BTC062900   | Shower Commode Seat - K Care - H/Made - Open 550mm - Blue                                       | \$ 579.00   | \$ 173.70   | 1   |
| BTC063100   | Shower Commode Seat - K Care - H/Made - Open 650mm  | \$ 1,029.00 | \$ 308.70   | 1   |
| BTC066300   | Bedpan and Bowl Carrier - 4HC (& KCare) - for Shower commodes                                   | \$ 210.00   | \$ 63.00    | 2   |
| BTC077000   | Pan Carrier - for K Care Junior 400mm Shower Commode - S/Steel                                  | \$ 293.00   | \$ 87.90    | 1   |
| BTC089124   | Shower Commode Seat - Aspire - Hand Made Pressure Reducing - 320mm - Closed Front               | \$ 346.00   | \$ 103.80   | 1   |
| BTS108604   | Transfer Disk Aquatec   | \$ 116.00   | \$ 34.80    | 1   |
| BTS151110   | Shower Chair - Saljol - Spa - Medium - Dark Grey - 150kg  | \$ 598.00   | \$ 179.40   | 1   |
| BTT156500   | Over Toilet Aid Toilet Seat Lid - for U Frame OTA   | \$ 23.70    | \$ 7.11     | 3   |
| BTT157120   | Throne Childrens Step   | \$ 402.00   | \$ 120.60   | 1   |
| BTT160600   | Urinal Uribag Female Foldaway Reusable  | \$ 57.00    | \$ 17.10    | 2   |
| CHP198975   | Chair Pressure Relief - Aspire COVE Chair - Navy/Bronx Steel - Large - 160kg                    | \$ 3,689.00 | \$ 1,844.50 | 1   |
| CHP213000   | Hilite Chair Tray Support Tube Vanilla  | \$ 265.00   | \$ 79.50    | 1   |
| CHP213300   | Hilite Chair Rear Adj Legs w/ Wheels - Vanilla/Champ  | \$ 201.00   | \$ 60.30    | 1   |
| CHP216355   | Chair Centre Directional Wheel - Air Comfort - Deluxe V2  | \$ 1,409.00 | \$ 422.70   | 1   |
| CHP227100   | Chair Powerlift Recline - Aspire Verona - Mushroom - Single Motor - 130kg ^                     | \$ 1,659.00 | \$ 829.50   | 1   |
| DIN265180   | Closure Cap 250ml (S.M.S.)  | \$ 2.00     | \$ 0.60     | 9   |
| DIN265780   | Tumbler Autoplas 140ml (S.M.S.)   | \$ 3.90     | \$ 1.17     | 54  |
| DIN266600   | Tumbler Autoplas 230ml Green  | \$ 5.49     | \$ 1.65     | 80  |
| DLG271880   | Plug Puller (S.M.S.)  | \$ 14.00    | \$ 4.20     | 4   |
| FBS281300   | Wrist Wrap - BodyAssist - Thermal Deluxe Wrap - Beige   | \$ 40.00    | \$ 12.00    | 1   |
| FBS283555   | Elastic Strap with Loop 70cm One Size White   | \$ 37.00    | \$ 11.10    | 1   |
| FBS283680   | Abdominal Support - Belt - BodyAssist - Small   | \$ 68.00    | \$ 20.40    | 1   |
| FBS284380   | Thigh Support - BodyAssist - Elastic Slip on - Sml  | \$ 57.00    | \$ 17.10    | 1   |
| FBS284670   | Elbow Sleeve - BodyAssist - Slip-On Elastic- Bge - XL   | \$ 24.90    | \$ 7.47     | 1   |
| FBS284990   | Wrist Strap - Adjustable - McDavid - Thermal - Black One-Size                                   | \$ 40.00    | \$ 12.00    | 1   |
| GMD333286   | Sphygmomanometer Bladder Latex Free Large Adult 1 Tube  | \$ 43.20    | \$ 12.96    | 1   |
| GMS336200   | Wipes White Rediwipe 300x330mm (Pkt 100)  | \$ 19.60    | \$ 5.88     | 2   |
| LSS371954   | Lifter Sling - Oxford - Access Loop Toileting - Polyester - Small                               | \$ 404.00   | \$ 202.00   | 1   |
| LSS375390   | Sling - Pivot Hygiene - Mesh Medium   | \$ 327.00   | \$ 163.50   | 1   |
| LSS380900   | Lifter Sling - Gen Purp Polyester - Kerry - Small   | \$ 689.00   | \$ 344.50   | 1   |
| LSS381000   | Lifter Sling - Gen Purp Polyester - Kerry - Medium  | \$ 609.00   | \$ 304.50   | 1   |
| LSS394150   | Sling - SupportVest w/ padding Mod 91 - Liko - XL - w/2 Loop Clips Fits Sabina II 350 SlingBar. | \$ 749.00   | \$ 224.70   | 1   |
| LSS394200   | Sling - ComfortVest Mod 95 - Liko - Lge - w/4 Loop Clips Fits Sabina II Comfort SlingBar.       | \$ 829.00   | \$ 248.70   | 1   |
| PED490550   | Seating Support - Firefly - GoTo Size 1 w/Advanced Headrest- Floorsitter- Cushion - Red Vinyl   | \$ 1,769.00 | \$ 1,238.30 | 1   |
| PTA493750   | HipSaver Slim Fit Veterans Kit Female X-Large   | \$ 529.00   | \$ 370.30   | 1   |
| PTA494001   | HipSaver Slim Fit Veterans Kit Female Extra Small   | \$ 529.00   | \$ 370.30   | 1   |
| PTA494150   | HipSaver Hip Protecting Pads (S.M.S.)   | \$ 64.68    | \$ 45.28    | 1   |
| PTA514100   | Abductor Button - Pride - Gel - Padded  | \$ 115.00   | \$ 80.50    | 1   |
| PTA516300   | Wedge Adductors - Jay - 1in - for Jay Cushions - Pair   | \$ 72.00    | \$ 50.40    | 1   |
| PTA533312L  | Limb Protection - DermaSaver - Zip Tube for Arms - Lge  | \$ 95.00    | \$ 66.50    | 1   |
| PTA533317XL | Limb Protection - DermaSaver - Leg Tube - XL  | \$ 115.00   | \$ 80.50    | 2   |
| PTA533328   | Limb Protection - Heel or Elbow - Foam - Sml  | \$ 79.00    | \$ 55.30    | 1   |
| PTC567800   | Cushion - JAY - Care 18x15.5in  | \$ 949.00   | \$ 474.50   | 1   |
| PTC587100   | Cushion - ROHO Nexus Spirit 18x16in   | \$ 726.00   | \$ 363.00   | 1   |
| PTC594920   | Cushion - ROHO Hybrid Elite 1R 16x18in Sensor Ready   | \$ 1,079.00 | \$ 539.50   | 1   |
| PTC601100   | Cushion - Synergy Simplicity 16 x 16in  | \$ 127.00   | \$ 63.50    | 1   |
| PTC607218   | Spex - Stnd-Cntr adj cushion 14in x 14in  | \$ 779.00   | \$ 389.50   | 1   |
| PTC621975   | Cushion - Axiom - Standard - 18 x 16  | \$ 157.00   | \$ 78.50    | 1   |
| PTM636625   | AirLo - Thermal Hyperstretch Microfibre Coverlet  | \$ 494.00   | \$ 148.20   | 1   |
| PWS655320   | Power Wheelchair - Whill - Model Ci - (Black Armrests) - 450mm - 114kg                          | \$ 7,489.00 | \$ 3,744.50 | 1   |
| TSE672526-S | Evacuation Sheet - Aidacare Lifecomfort - Single - 1900mm x 720mm - 300kg                       | \$ 150.00   | \$ 75.00    | 1   |

# WAREHOUSE CLEARANCE

ALL PRODUCTS  
BETWEEN  
**50-70%**  
OFF

HURRY,  
WHILST  
STOCK  
LASTS!

| Item Number | Item Description   | WAS       | NOW       | QTY |
|-------------|--|-----------|-----------|-----|
| WAA686250   | Crutch Acc - Aspire - Grey Padded Hand Grip                          | \$ 20.70  | \$ 6.21   | 1   |
| WAA687000   | Walker Acc - Basket for Seat Walker                                  | \$ 31.00  | \$ 9.30   | 1   |
| WAF705600BL | Seat Walker - Aspire PUSH DOWN - Blue - 7in Wheels - 130kg           | \$ 207.00 | \$ 103.50 | 1   |
| WAF705600RE | Seat Walker - Aspire PUSH DOWN - Red - 7in Wheels - 130kg            | \$ 207.00 | \$ 103.50 | 1   |
| WAF750150   | Seat Walker - Freedom - Stride - Red - 135kg                         | \$ 369.00 | \$ 184.50 | 1   |
| WAF751100   | Seat Walker - Saljol - Page Indoor Walker - Small - Charcoal - 150kg | \$ 380.00 | \$ 114.00 | 2   |
| WAS748110   | Walking Stick - Freedom - Cobra Cane - 710-940mm - 110kg             | \$ 56.00  | \$ 16.80  | 1   |
| WAS748140   | Walking Stick - Freedom - Crook Handle - Black - 770-970mm - 110kg   | \$ 28.50  | \$ 8.55   | 1   |
| WAS748160   | Walking Stick - Freedom - Crook Handle - Bronze - 770-970mm - 110kg  | \$ 28.50  | \$ 8.55   | 1   |
| WAS748180   | Walking Stick - Freedom - Crook Handle - Floral - 770-970mm - 110kg  | \$ 28.50  | \$ 8.55   | 4   |
| WAS748250   | Walking Stick - Freedom - Offset Handle - Black - 770-970mm - 110kg  | \$ 24.00  | \$ 7.20   | 4   |
| WAS748260   | Walking Stick - Freedom - Offset Handle - Bronze - 770-970mm - 110kg | \$ 24.00  | \$ 7.20   | 6   |
| WAS748270   | Walking Stick - Freedom - Offset Handle - Chrome - 770-970mm - 110kg | \$ 24.00  | \$ 7.20   | 4   |
| WAS748280   | Walking Stick - Freedom - Offset Handle - Tartan - 770-970mm - 110kg | \$ 24.00  | \$ 7.20   | 4   |
| WAS748380   | Walking Stick - Freedom - Swan Neck - Chrome - 770-970mm - 110kg     | \$ 28.50  | \$ 8.55   | 2   |
| WAS748390   | Walking Stick - Freedom - Crook Handle - Chrome - 770-970mm - 110kg  | \$ 28.50  | \$ 8.55   | 1   |
| WAS748450   | Walking Stick - Freedom - T Handle Ergo - Chrome - 770-970mm - 110kg | \$ 25.20  | \$ 7.56   | 5   |
| WAS748460   | Walking Stick - Freedom - T Handle Ergo - Leaves - 770-970mm - 110kg | \$ 24.00  | \$ 7.20   | 2   |