

WAREHOUSE CLEARANCE

ALL PRODUCTS BETWEEN **50-70% OFF**

HURRY, WHILST STOCK LASTS!

Item Number	Item Description	WAS	NOW	QTY
BEB045850	Bed - Aidacare AC3 V2 - Single - 3 Function - Teak Bed Ends	\$ 2,450.00	\$ 1,225.00	4
BEB045567	Bed - Invacare - CS7-42" Wide Deck (bed ends sep)	\$ 5,659.00	\$ 2,829.50	4
MDC409000	Catheter Y-Suction 12Fg Mini Pack - Pkt/50(old Code 110.12.580)	\$ 179.00	\$ 53.70	200
LIC350240	Ceiling Hoist - Prism Transactive TXE270D6B (C625) - Pwr Traverse w/6 Button Rtrn to Charge - 270kg	\$ 8,979.00	\$ 2,693.70	1
CHP199500	Chair Powerlift Recline - Air Comfort - Airwing - Lge 600x510x710mm - 1 Motor	\$ 3,829.00	\$ 1,148.70	3
CHP199600	Chair Powerlift Recline - Air Comfort - Airwing - Med 500x510x710mm - 1 Motor	\$ 3,829.00	\$ 1,148.70	2
CHP198970	Chair Pressure Relief - Aspire COVE Chair - Grey/Andreas Pumice - Large - 160kg ^	\$ 3,689.00	\$ 1,844.50	2
CHP198965	Chair Pressure Relief - Aspire COVE Chair - Navy/Bronx Steel - Medium - 160kg ^	\$ 3,689.00	\$ 1,844.50	2
WAT905150	Emerg Cart TRIBUTE 110 PKG - 10H DR321 - BRK AWY Lock & RHS Handle - CPR Bd, IV, CO2, Defib, Waste -	\$ 3,829.00	\$ 1,914.50	3
LIC360760	FreeSpan Rail 5m - with Likorall Trolley - 200kg	\$ 3,049.00	\$ 1,524.50	2
BEB050201	Hospital Bed - Malvestio Botero Bariatric Ward Bed w/Integrated Weigh Scale - 400kg	\$ 21,989.00	\$ 10,994.50	1
TSE674210	Hospital Stretcher - Aspire Advanced Emergency Stretcher - Hydraulic - 2140x840x590-900mm - 225kg	\$ 8,309.00	\$ 4,154.50	1
WAT905440	Paed Cart TRIBUTE 110 PKG - 10H DR810 - BRK AWY Lock & RHS Handle - CPR Bd, IV, CO2, Defib, Waste (\$ 4,179.00	\$ 2,089.50	2
PWC655350	Power Wheelchair - Quickie Q500 - Front Wheel Drive - (as per script)	\$ 17,819.00	\$ 8,909.50	1
CHT100210	Treatment Chair - Greiner Multiline Next DC - Navy Blue Vinyl - 220kg	\$ 8,789.00	\$ 2,636.70	1
CHT100100	Treatment Chair - Greiner Multiline Next IT - Navy Blue Vinyl - 220kg	\$ 6,500.00	\$ 1,950.00	1
CHT100110	Treatment Chair - Greiner Multiline Next IT V2 - Navy Blue Vinyl - 220kg	\$ 8,249.00	\$ 2,474.70	2
LIC360860	UltraTwin FreeSpan Upright Support - Incl. wheels, charger plate and parking bracket for sling bar	\$ 6,399.00	\$ 1,919.70	1
LSS363095	Voyager Portable Hoist	\$ 5,599.00	\$ 1,679.70	1
MWC427331	Wheelchair - Quickie Xenon 2 (As per script)	\$ 4,989.00	\$ 2,494.50	1

Item Number	Item Description	WAS	NOW	QTY
BEA007488	Alarm - Cura - Cordless Floor Mat - (1220 x 610) includes CU-2619 Transmitter	\$ 649.00	\$ 194.70	1
BEA007850	Alarm - Sofihub TEQ Home (1 x SOFIHUB TEQ Home + 8 Sensors)	\$ 1,419.00	\$ 709.50	1
BEA013800	Bed Pole - Dual Prong for Single Bed	\$ 52.00	\$ 15.60	1
BEA014605	Bed Cradle - Aspire - Fixed - Treated Steel ^	\$ 69.00	\$ 20.70	8
BEA024100	Bed End Set - CS7 - Amherst 36in w/ Ctrl Panel - Solar Oak	\$ 579.00	\$ 173.70	1
BEA028100	Overbed/Chair Table - Aspire - Wide C Base - Thermoform Recessed Split Tilt Top - Cream	\$ 349.00	\$ 174.50	2
BEA029845	Bed End Shrouds - Beech - Single - Aspire COMMUNITY CARE Bed (Premium)	\$ 300.00	\$ 90.00	4
BEA824130	Bed Moving Trolley - Generic	\$ 1,029.00	\$ 308.70	2
BEA826110	Folding Side Rails - 3/4 length - Aspire Elite Ward Bed	\$ 199.00	\$ 59.70	2
BEB045567	Bed - Invacare - CS7-42" Wide Deck (bed ends sep)	\$ 5,659.00	\$ 2,829.50	4
BEB045850	Bed - Aidacare AC3 V2 - Single - 3 Function - Teak Bed Ends	\$ 2,450.00	\$ 1,225.00	4
BEB045860	Bed - Aidacare AC3 V2 - King Single - 3 Function - Teak Bed Ends	\$ 2,950.00	\$ 1,475.00	1
BEB050201	Hospital Bed - Malvestio Botero Bariatric Ward Bed w/Integrated Weigh Scale - 400kg	\$ 21,989.00	\$ 10,994.50	1
BEM046450	Mattress - Aspire Lifecomfort - Standard - Narrow - 2000x850x130mm - 120kg	\$ 519.00	\$ 259.50	1
BEM046500	Mattress - Aspire Lifecomfort - Deluxe - Single - 1940x900x150mm - 150kg	\$ 589.00	\$ 294.50	1
BEM046725	Mattress - Aspire Lifecomfort Premier Care - Single - 180kg	\$ 869.00	\$ 434.50	2
BEP051620	Pillow - Travel Nut - Charcoal	\$ 40.20	\$ 20.10	1
BTC061800	Shower Commode - K Care - Platform Footrest - Maxi 600mm 210kg - S/Steel (seat sep)	\$ 3,209.00	\$ 1,604.50	1
BTC062400	Shower Commode Seat - K Care - H/Made - Closed 500mm	\$ 749.00	\$ 224.70	3
BTC062500	Shower Commode Seat - K Care - H/Made - Closed 550mm	\$ 839.00	\$ 251.70	1
BTC066300	Bedpan and Bowl Carrier - 4HC (& KCare) - for Shower commodes	\$ 210.00	\$ 63.00	1
BTC066500	Shower Commode Seat - Aspire - Closed 410mm	\$ 143.00	\$ 42.90	1
BTC067108S	Shower Commode Elevating Legrest - Aspire - LEFT - Non-Tilt Commodes Only - Shiny Finish	\$ 370.00	\$ 111.00	1
BTC067109S	Shower Commode Elevating Legrest - Aspire - RIGHT - Non-Tilt Commodes Only - Shiny Finish	\$ 370.00	\$ 111.00	1
BTC071000	Shower Commode - K Care - SA Footrest - 400mm - S/Steel (seat sep)	\$ 2,079.00	\$ 623.70	1
BTC074000	Shower Commode Seat - K Care - H/Made - Open Left Side	\$ 899.00	\$ 269.70	1
BTC075300	Backrest - K Care - Curved Padded - STD Width - for K Care Shower Commodes - Fitted	\$ 310.80	\$ 93.24	1
BTC076600	Bowl Carrier - for K Care Shower Commodes - S/Steel	\$ 168.00	\$ 50.40	1
BTC078300	Anti Tip Bars - for K Care Shower Commodes - OLD - Compatible w Candy Cane Leg	\$ 165.60	\$ 49.68	1
BTC085200	Shower Commode - K Care - Child - Self Prop - 320mm (seat sep)	\$ 3,119.00	\$ 935.70	1
BTC087640	Commode Seat - Invacare - Universal Soft - for Aquatec Ocean Shower Commodes	\$ 237.00	\$ 71.10	1
BTC087648	Insert - Invacare - for Aquatec Shower Commodes - Recess Cover - Soft Seat	\$ 116.00	\$ 34.80	1
BTC089191	Shower Commode Seat - Aspire Custom 4HC - Hand Made Pressure Reducing - 550mm - Open Front	\$ 461.00	\$ 138.30	1
BTR101420	Multi Angle Fittings for Tubes	\$ 47.10	\$ 14.13	6
BTR101751	Rail 1500mm Tube 32mm diameter BGE	\$ 73.80	\$ 22.14	11
BTR101771	One T-piece, one wall flange with washer & screw BGE	\$ 19.10	\$ 5.73	8
BTR101790	One 45 deg bend, one wall flange with fixings	\$ 31.40	\$ 9.42	5
BTR101801	One 90 deg bend, one wall flange with fixings BGE	\$ 31.40	\$ 9.42	3

WAREHOUSE CLEARANCE

ALL PRODUCTS BETWEEN **50-70% OFF**

HURRY, WHILST STOCK LASTS!

Item Number	Item Description	WAS	NOW	QTY
BTR101830	4 way fitting, wall flange includes fixings	\$ 31.10	\$ 9.33	2
BTR101883	Plug for End of Tube	\$ 3.00	\$ 0.90	2
BTR101885	Internal Corner fitting, one corner wall flange BGE	\$ 35.50	\$ 10.65	7
BTR101967	90 degree tube connector not elbow	\$ 16.40	\$ 4.92	6
BTR102930	Shower Slide Rail	\$ 31.40	\$ 9.42	1
BTS106000	Shower Stool w/ Arms & Padded Seat - K Care - 550mm Seat - Zinc	\$ 778.40	\$ 233.52	3
BTS106900	Bath Chair Acc - Myco Swivel Bather Spacer 6in (150mm)	\$ 142.00	\$ 42.60	1
BTS112253	Bath Lift Accessory - Aquatec Orca - Pommel 7cm diameter	\$ 277.00	\$ 83.10	1
BTS113600	Shower Stool w/ Arms - Dayz - PVC Strapping Seat - Zinc	\$ 145.00	\$ 43.50	1
BTS116100	Shower Chair - Wall Mounted - Fold Away - K Care - AS Large Padded 125kg	\$ 2,358.57	\$ 707.57	1
BTS119300	Shower Stool - in U Frame Adjustable	\$ 97.00	\$ 29.10	1
BTT149700	Toilet Support Rail Toilet Roll Holder - for Devon Rail	\$ 74.00	\$ 22.20	1
BTT156700	Over Toilet Aid - U Frame Adjustable - 510mm (splashguard sep)	\$ 112.90	\$ 33.87	1
CHP191200	Chair Highback - K Care Throne Chair - 700mm - Fixed Height - 400kg	\$ 3,229.00	\$ 1,614.50	1
CHP191510	Chair Lowback - K Care Utility - 650mm - Grey/Slate	\$ 2,189.00	\$ 1,094.50	1
CHP197055	Chair Powerlift Recline - Pride - C101 Riser Recliner Pecan - 510x510x700mm - 1 Motor	\$ 1,391.50	\$ 695.75	1
CHP197741	Chair Powerlift Recline - Aspire Signature - Dual Action - Small - Heritage Vinyl - 155kg	\$ 2,269.00	\$ 1,134.50	2
CHP197770	Head/Armrest Cover Set - Aspire Signature Recliners - Mink	\$ 46.30	\$ 13.89	2
CHP197772	Head/Armrest Cover Set - Aspire Signature Recliners - Heritage	\$ 46.30	\$ 13.89	2
CHP197825	Variable Angle Leg Rest - Atlanta & Monaco Chairs	\$ 519.00	\$ 155.70	1
CHP197865	Leg Laterals - Phoenix & Sorrento Chairs	\$ 639.00	\$ 191.70	1
CHP197873	Lateral Supports - Seating Matters - Denver Chairs	\$ 679.00	\$ 203.70	2
CHP198810	Overnight Chair/Bed - Aspire SLEEPER Chair - Beige Vinyl - 150kg	\$ 3,109.00	\$ 1,554.50	1
CHP198820	Overnight Recliner Chair - Aspire RECOVERY Chair - Beige Vinyl - 150kg	\$ 1,109.00	\$ 554.50	1
CHP198960	Chair Pressure Relief - Aspire COVE Chair - Grey/Andreas Pumice - Medium - 160kg ^	\$ 3,689.00	\$ 1,844.50	1
CHP198965	Chair Pressure Relief - Aspire COVE Chair - Navy/Bronx Steel - Medium - 160kg ^	\$ 3,689.00	\$ 1,844.50	2
CHP198970	Chair Pressure Relief - Aspire COVE Chair - Grey/Andreas Pumice - Large - 160kg ^	\$ 3,689.00	\$ 1,844.50	2
CHP198975	Chair Pressure Relief - Aspire COVE Chair - Navy/Bronx Steel - Large - 160kg	\$ 3,689.00	\$ 1,844.50	1
CHP199500	Chair Powerlift Recline - Air Comfort - Airwing - Lge 600x510x710mm - 1 Motor	\$ 3,829.00	\$ 1,148.70	3
CHP199600	Chair Powerlift Recline - Air Comfort - Airwing - Med 500x510x710mm - 1 Motor	\$ 3,829.00	\$ 1,148.70	2
CHP204500	Chair Highback - K Care Hilite - 490mm - Grey/Slate	\$ 560.00	\$ 168.00	1
CHP208500	Chair Lowback - K Care Utility - 440mm - Vanilla/Champ	\$ 355.00	\$ 106.50	1
CHP210700	Chair Leg Rest - Adjustable - Champagne Vinyl	\$ 292.00	\$ 87.60	2
CHP213200	Hilite Chair Rear Adj Legs w/ Wheels - Grey/Slate	\$ 201.00	\$ 60.30	1
CHP213400	Stool - Kitchen Stool - Adjustable No Arms - 360mm	\$ 509.00	\$ 152.70	3
CHP216100	Chair Pressure Relief - Air Comfort - Princess Bed - Boston	\$ 3,039.10	\$ 911.73	1
CHP216300	Chair Footrest - Princess Chair/Bed & Deluxe - Air Comfort	\$ 477.00	\$ 143.10	1
CHP216352	Chair Footrest - Air Comfort - Deluxe V2 - Large old style	\$ 331.00	\$ 99.30	3
CHP216400	Chair Tray Small - Princess Chair/Bed - Air Comfort	\$ 149.00	\$ 44.70	1
CHP217300	Chair Tray Large - Tilt / Deluxe Bed/Daily Chair - Air Comfort	\$ 262.00	\$ 78.60	1
CHP227320	Chair Powerlift Recline - Aspire Raphael Quattro - MAXI - Bronze - 4 Motor - 220kg	\$ 3,149.00	\$ 944.70	1
CHT100100	Treatment Chair - Greiner Multiline Next IT - Navy Blue Vinyl - 220kg	\$ 6,500.00	\$ 1,950.00	1
CHT100110	Treatment Chair - Greiner Multiline Next IT V2 - Navy Blue Vinyl - 220kg	\$ 8,249.00	\$ 2,474.70	2
CHT100210	Treatment Chair - Greiner Multiline Next DC - Navy Blue Vinyl - 220kg	\$ 8,789.00	\$ 2,636.70	1
CHT100400	Treatment Chair - Aspire DIALYSIS Chair - Beige Vinyl - 150kg	\$ 1,869.00	\$ 560.70	2
CON250000	TENA Pads for Men Level 2 (Pkt 10)	\$ 19.90	\$ 5.97	2
DCO100225	Dementia Care -Home Safety - Cupboard C-Lock	\$ 24.90	\$ 7.47	1
DLK276400	Cutlery - Amefa - Contoured Fork Angled Right	\$ 52.31	\$ 15.69	1
DLK276410	Cutlery - Amefa - Contoured Spoon Angled Left	\$ 93.00	\$ 27.90	1
DLK278110	Neater Eater Carry Bag - for Robotic Neater Eater	\$ 364.00	\$ 109.20	1
FBS280020	Actimove Sling Collar & Cuff Support (2 x Rolls 12m)	\$ 168.00	\$ 50.40	1
FBS280370	Abdominal Binder - Orthobinder - 6in - Extra Small - 44-58cm (17.5-23in)	\$ 82.00	\$ 24.60	4
FBS281202	Foot Support - BodyAssist - Innersoles - Silicone Gel - Full Length - Med	\$ 56.00	\$ 16.80	1
FBS285450	Back Support Obus Forme High	\$ 235.00	\$ 70.50	1
FBS291350	Splint Pan Netting	\$ 52.00	\$ 15.60	1
FBS291450	Norco Thumb Wrap with Wrist Support Right	\$ 75.00	\$ 22.50	1
FBS301170S2	Hand - Gloves - Handeze Support Gloves - Pair - Med - Bge	\$ 63.00	\$ 18.90	1
FBS301258	Shoulder Support - Rehband - EQ posture support - X Sml	\$ 131.00	\$ 39.30	1
FBS301269M	Abdominal Binder - Oppo - Firm - Med	\$ 62.00	\$ 18.60	1
FBS301279	Lumbar - Oppo - Foam Sacro Lumbar Suppt - Sml	\$ 87.00	\$ 26.10	1
FBS301520	Knee Brace - Rehband - Genu Carezza Patella Stabilizer Wrap - O/Patella - XX Sml	\$ 62.80	\$ 18.84	1
FBS3015803B	Shoe - Post Op - Darco - Orthowedge Shoe - Lge	\$ 97.00	\$ 29.10	1

WAREHOUSE CLEARANCE

ALL PRODUCTS BETWEEN **50-70% OFF**

HURRY, WHILST STOCK LASTS!

Item Number	Item Description	WAS	NOW	QTY
GMD334050	Oxygen Watch - Pulse Oximeter	\$ 248.00	\$ 74.40	1
GSS340340	BetterLiving Cushion Wheelchair Lumbar 20"	\$ 86.00	\$ 25.80	1
LAC353030	Linen Bag Waterproof Black	\$ 70.00	\$ 21.00	2
LIC350240	Ceiling Hoist - Prism Transactive TXE270D6B (C625) - Pwr Traverse w/6 Button Rtrn to Charge - 270kg	\$ 8,979.00	\$ 2,693.70	1
LIC353100	Ceiling Hoist - Prism - Track - 2.0m Rail	\$ 176.00	\$ 52.80	1
LIC353235	Ceiling Hoist - Prism - Track w/Wings - 400mm - Maintenance Openinig	\$ 125.00	\$ 37.50	38
LIC353550	Ceiling Hoist - Prism - Trolley for Portable Hoists - Fixed Eyelet	\$ 101.00	\$ 30.30	2
LIC354115	Soft Strap - P200 Weigh Scale	\$ 58.00	\$ 17.40	1
LIC355871	Back Cover Only - Prism - For P200	\$ 37.20	\$ 11.16	1
LIC357057	Ceiling Hoist Connector Bracket 150mm	\$ 48.60	\$ 14.58	16
LIC357310	Track with Wings 3.0m	\$ 393.00	\$ 117.90	2
LIC359410	In Rail Charge Carriage - LikoGuard - (with Step up card and cable) - Max 363kg (800lbs)	\$ 559.00	\$ 167.70	1
LIC359420	In-Rail-Charge Charger - LikoGuard, AU	\$ 462.00	\$ 138.60	1
LIC359490	Ceiling Hoist - Likorall 200 - natural - incl hand control - 200 kg (440 lbs) SWL	\$ 2,929.00	\$ 878.70	1
LIC359560	Hand Control - Likorall ES-4MT (for Lifting and Transfer Motor)	\$ 483.00	\$ 144.90	3
LIC359570	Hand Control - Likorall ES-6MTS (for Lifting, Transfer Motor and Switches)	\$ 579.00	\$ 173.70	2
LIC360320	Pendant Cover dia 50, 2-pack (2.0 in.) - Liko	\$ 13.80	\$ 4.14	16
LIC360630	Side Rail Switch, electrical - White - Liko	\$ 1,749.00	\$ 524.70	1
LIC360670	HandControl HB80 - Liko	\$ 374.00	\$ 112.20	2
LIC360680	Extension cable HB80 - Liko	\$ 287.00	\$ 86.10	5
LIC360720	Traverse Carrier Multi, wide - Incl. 1 pce Traverse Rail Carrier, wide. Fits Likorall ES.	\$ 1,329.00	\$ 398.70	2
LIC360760	FreeSpan Rail 5m - with Likorall Trolley - 200kg	\$ 3,049.00	\$ 1,524.50	2
LIC360860	UltraTwin FreeSpan Upright Support - Incl. wheels, charger plate and parking bracket for sling bar	\$ 6,399.00	\$ 1,919.70	1
LIC360900	UltraTwistTM Wide - incl, hand control. Max 500 kg (1100 lbs).	\$ 3,829.00	\$ 1,148.70	1
LIC360940	Pendant Plate 90 - E-System - Liko	\$ 124.00	\$ 37.20	4
LIC360980	Curtain Rail T, 2.5m - T-System - Liko	\$ 569.00	\$ 170.70	2
LSS363095	Voyager Portable Hoist	\$ 5,599.00	\$ 1,679.70	1
LSS367750	Sling Liko Original H/B Thin Net Poly Reinf Leg Support w Belt	\$ 442.00	\$ 221.00	1
LSS368400	Sling Liko Hygiene Nylon Reinforced Leg Support Large	\$ 349.00	\$ 174.50	1
LSS370100	Leg Harness Liko Medium/Large	\$ 224.00	\$ 112.00	1
LSS371642	Quickfit Deluxe H/S Pad Large	\$ 489.00	\$ 244.50	1
LSS371699	Lifter Sling - Oxford - Quickfit Deluxe Poly - Large	\$ 283.00	\$ 141.50	1
LSS374905	Sling - Pivot - Mesh - XXL	\$ 261.00	\$ 130.50	1
LSS377763	Lifter Sling - Allegro - Head Support Pivot Fabric - Large	\$ 444.00	\$ 222.00	1
LSS377774	Lifter Sling - Allegro - Standing Keyhole - Fabric - XLarge	\$ 619.00	\$ 309.50	1
LSS381000	Lifter Sling - Gen Purp Polyester - Kerry - Medium	\$ 609.00	\$ 304.50	1
LSS382400	Lifter Sling - Full Body Polyester w/ Div Leg - Kerry - Medium LR408N2	\$ 376.00	\$ 188.00	1
LSS386202	Lifter - Invacare - 2 Point Spreader Bar - Jasmine	\$ 262.00	\$ 78.60	1
LSS388500	Transfer - Board - Patient Slide Large 1500x580mm	\$ 589.00	\$ 176.70	1
LSS390960	Lifter Sling - Aspire - SML - Low Friction General Purpose w/Head Supp - CLIPS - Parasilk - 200kg	\$ 539.00	\$ 269.50	1
LSS393410	Wall Mounted Battery Charger, AU/NZ for CBL Battery.	\$ 679.00	\$ 203.70	1
LSS393540	Spreader Bar - Liko Universal SlingBar 600 - Max 300 kg (660 lbs).	\$ 490.00	\$ 147.00	1
LSS393610	Sling Cross-bar 670 with Quick-release Hook - Liko - Max 300kg	\$ 1,019.00	\$ 305.70	1
LSS393800	Adapter 22 for Quick-release Hook - Liko - Fits Liko MasterLift from -89 and Likorall.	\$ 53.00	\$ 15.90	1
LSS394120	Sling - SupportVest w/ padding Mod 91 - Liko - Sml - w/2 Loop Clips Fits Sabina II 350 SlingBar.	\$ 709.00	\$ 212.70	1
LSS394170	Sling - SafetyVest Mod 93 - Liko - Lge - w/4 Loop Clips Fits Sabina II 350 SlingBar	\$ 849.00	\$ 254.70	1
LSS394190	Sling - ComfortVest Mod 95 - Liko - Med - w/4 Loop Clips Fits Sabina II Comfort SlingBar.	\$ 799.00	\$ 239.70	1
LSS394210	Sling - ComfortVest Mod 95 - Liko - XL - w/4 Loop Clips Fits Sabina II Comfort SlingBar.	\$ 849.00	\$ 254.70	1
LSS394220	Sling - HygieneVest Highback Mod 55 - Liko - Poly w/ reinforced leg supports - Sml - 200kg	\$ 1,029.00	\$ 308.70	1
LSS394709	Sling - Liko - Original MOD 10/11 - Net Polyester -M	\$ 284.00	\$ 142.00	1
MDC409000	Catheter Y-Suction 12Fg Mini Pack - Pkt/50(old Code 110.12.580)	\$ 179.00	\$ 53.70	200
MDC414860	Bag ULTIMATE "T" BOX/100	\$ 392.00	\$ 117.60	5
MWC427331	Wheelchair - Quickie Xenon 2 (As per script)	\$ 4,989.00	\$ 2,494.50	1
MWC449214	Wheelchair - Start M2S - Standard - 405mm Wide (as per script)	\$ 949.00	\$ 474.50	1
MWS439601	WChair Folding Std - Karma Ergo Lite Deluxe - Self Prop	\$ 1,239.00	\$ 619.50	1
MWS442370	Wheelchair Folding - Aspire METRO X - Self Prop - METALLIC ORANGE - 100kg ^	\$ 679.00	\$ 339.50	2
MWS442800	Wheelchair Folding Push/Pull - 500mm Leverdrive - Left Hand	\$ 4,629.00	\$ 2,314.50	1
MWS449600	WChair - Tilt in Space - Aspire Rehab RS - 410mm - 24in Wheels - 180kg - Gunmetal Grey	\$ 2,479.00	\$ 1,239.50	1
MWS450200	Ward Transport Chair - Aspire Swish HD - 740mm - 300kg	\$ 1,995.00	\$ 997.50	1
PED478453	Kimba Headrest Supports Size 1	\$ 559.00	\$ 391.30	1
PED480807	Paediatric Walking Aid - Nimbo Lightweight - Adult 920-1040mm	\$ 529.00	\$ 370.30	1
PED480820	Paediatric Walking Aid - Trekker Gait Trainer (as per script)	\$ 769.00	\$ 538.30	1

WAREHOUSE CLEARANCE

ALL PRODUCTS BETWEEN **50-70% OFF**

HURRY, WHILST STOCK LASTS!

Item Number	Item Description	WAS	NOW	QTY
PED480829	Paediatric Bathing Aid - Otter Bathing System - Small	\$ 779.00	\$ 545.30	1
PED480965	Seat Harness - Numbo - for Tyke to Youth 390 to 770mm	\$ 205.00	\$ 143.50	1
PED490220	Hospital Bed - Aspire Paediatric Cot & Bath - Small - 700x500x610-860mm	\$ 695.00	\$ 486.50	1
PTA508615	WChair Footplate Protector Cover - Pride - 27.94 x 22.86cm - Gel - Single Piece	\$ 231.00	\$ 161.70	1
PTA508810	WChair Knee Protector Buttons - Pride - 1.98cm - for Swing Away Legrests - Gel	\$ 203.00	\$ 142.10	1
PTA511000	HipSaver - Nursing Home - Starter Kit (3-Pants, 1-Pads) - Male - Sml	\$ 269.00	\$ 188.30	1
PTA513914	HipSaver - QuickChange - Starter Kit (3-Pants, 1-Pads) - Male - Med	\$ 269.00	\$ 188.30	1
PTA521704	Boot - Sheepskin Hi Temp - Wrap Around - C/Toe - XX Sml	\$ 191.00	\$ 133.70	1
PTA521898	Boot - Sheepskin Hi Temp - Short Slipper - C/Toe - X Sml	\$ 215.00	\$ 150.50	2
PTA521904	Sheepskin Hi Temp Sovereign Snug Extra Small (S.M.S.)	\$ 121.35	\$ 84.95	3
PTA532965	DermaSaver Shin Tube Medium	\$ 77.00	\$ 53.90	1
PTA533317XL	Limb Protection - DermaSaver - Leg Tube - XL	\$ 115.00	\$ 80.50	4
PTA533318XL	Limb Protection - DermaSaver - Leg Tube with Double Knee - XL	\$ 132.00	\$ 92.40	1
PTA534101	Limb/Cast Protection - Limbo - Adult Hand - Waterproof	\$ 51.00	\$ 35.70	1
PTB532360	WChair Back - JAY J2 - 18in Tall	\$ 1,129.00	\$ 564.50	1
PTB533705	WChair Back - JAY J3 - Posterior - Lower Thoracic Med 18in	\$ 1,069.00	\$ 534.50	1
PTB533838	WChair Back - JAY J3 - Laterals - Swing Away 6x6in	\$ 759.00	\$ 379.50	1
PTB535064	Vigour-Mid Active Back Support for 16in W/C (16in high) for mid trunk support	\$ 1,329.00	\$ 664.50	1
PTB535406	SuperShape Back Support - 16in Wide x 18in High with adjustable shape and contour system	\$ 2,609.00	\$ 1,304.50	1
PTB539870	WChair Back - Dreamline Support S3 Backrest - Quick Release - 18in x 22in	\$ 1,529.00	\$ 764.50	1
PTC536149	J3 Field Adjustable Supplement Fluid Pad Small	\$ 64.00	\$ 19.20	1
PTC536150	J3 Field Adjustable Supplement Fluid Pad Large	\$ 73.00	\$ 21.90	1
PTC551285	Cushion - ROHO AirLite 14 x 14in	\$ 389.00	\$ 194.50	1
PTC551580	Cushion Vicair Adjuter 10 & incotec Cover - 16x16in	\$ 899.00	\$ 449.50	3
PTC555100	Cushion - ROHO Harmony 20x18in	\$ 409.00	\$ 204.50	1
PTC555200	Cushion Cover - ROHO Harmony 16x16in	\$ 143.00	\$ 42.90	1
PTC556600	Cushion Cover - ROHO Waterproof HP 16.75x18.5in	\$ 179.00	\$ 53.70	3
PTC561274	Cushion - Jay X2 - 17 x17in	\$ 759.00	\$ 379.50	4
PTC568000	Cushion - JAY - Care 16x15.5in	\$ 949.00	\$ 474.50	1
PTC571913	Cushion - JAY - Lite 16x18in	\$ 509.00	\$ 254.50	1
PTC572800	Cushion - JAY - Combi Soft 20x18in	\$ 367.00	\$ 183.50	1
PTC576100	Cushion Cover - JAY J2 - 20x20in - Recline Incontinence	\$ 186.00	\$ 55.80	1
PTC580050	Cushion - JAY - Easy 15x15in Curved Bottom	\$ 779.00	\$ 389.50	2
PTC585800	Cushion - ROHO Mini-Max 15x15in	\$ 959.00	\$ 479.50	1
PTC586700	Cushion - ROHO Nexus Spirit 15x15in	\$ 726.00	\$ 363.00	1
PTC594920SR	Cushion - ROHO Hybrid Elite 1R 16x18in w/HD cover Sensor Ready	\$ 1,079.00	\$ 539.50	1
PTC594929	Cushion - ROHO Hybrid Elite 1R 17x16in Sensor Ready	\$ 1,079.00	\$ 539.50	1
PTC594954	Cushion - ROHO Hybrid Elite 2R 18x18in	\$ 1,079.00	\$ 539.50	1
PTC595130	Cushion - ROHO Hybrid Elite Cover 18x20in	\$ 197.00	\$ 59.10	1
PTC595200	ROHO Cover 8x9 LP (S.M.S.)	\$ 89.00	\$ 26.70	1
PTC595800	Cushion Cover - ROHO Nexus Spirit 16x18in	\$ 179.00	\$ 53.70	1
PTC596150	Cushion Cover - ROHO Universal HP 18.5x16in	\$ 153.00	\$ 45.90	1
PTC598300	Cushion Invacare Ulitmate 18x18"	\$ 1,069.00	\$ 534.50	1
PTC598700	Cushion - Synergy Solution 16 x 16in	\$ 199.00	\$ 99.50	1
PTC607092	Spex - Kidz High-Contour adj cushion 13in x cstm w/	\$ 879.00	\$ 439.50	1
PTC607115	Vigour - Kidz Standard Contour Cushion 13inx14in	\$ 649.00	\$ 324.50	1
PTC607174	Vigour - Standard Contour Cushion 16inx18in	\$ 749.00	\$ 374.50	1
PTC607311	Spex Super High-Contour adjustable cushion 14in x 14in	\$ 1,099.00	\$ 549.50	1
PTC620610	Cushion - Dreamline Support Cushion - 16in x 20in	\$ 929.00	\$ 464.50	1
PTC621200	Cushion - Dreamline G3 Cushion - Foam - 18in x 22in	\$ 1,259.00	\$ 629.50	1
PTC622130	Cushion - Pride Stealth Solution Cushion - 12in x 12in	\$ 579.00	\$ 289.50	1
PTC622150	Cushion - Pride Stealth Solution Cushion - 14in x 14in	\$ 579.00	\$ 289.50	1
PTC622270	Cushion - Pride Stealth Solution Cushion - 17in x 20in	\$ 579.00	\$ 289.50	1
PTC622310	Cushion - Pride Stealth Solution Cushion - 18in x 20in	\$ 579.00	\$ 289.50	1
PTC622430	Cushion - Pride Stealth Solution Cushion - 22in x 22in	\$ 669.00	\$ 334.50	1
PTC622790	Cushion - Pride Stealth Spectrum Cushion - Foam - 26in x 24in	\$ 749.00	\$ 374.50	1
PTC623400	Cushion - Pride Stealth Spectrum Cushion - Gel - 22in x 22in	\$ 739.00	\$ 369.50	1
PTC623600	Cushion - Pride Stealth Tru Comfort Cushion - 20in x 20in	\$ 629.00	\$ 314.50	1